
DK	Ståsejl og Hoftestøttesejl	2
GB/US	Standing Sling and Hip Support Sling	6
DE	Standstz und Hüftstütze	10
SE	Ståsele och stödsele för höften	14
NO	Ståseil og hofteseil	18
FR	Harnais de station debout et harnais de maintien des hanches	22
IT	Imbragatura per posizione eretta e imbragatura di supporto dell'anca	26
NL	Corsetband en Heupband	30

User manual – vers. 103.0

Guldmann™

DK... STÅSEJL OG HOFTESTØTTESEJL

Vers. 103.0

Varenr.:

Ståsejl

28336 S/M

28337 L

28339 XL

28340 XXL

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28338 One size

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GB... STANDING SLING AND HIP SUPPORT SLING

Vers. 103.0

Item nos:
 Standing Sling
 28336 S/M
 28337 L
 28339 XL
 28340 XXL

Hip Support Sling
 28338 One size

1.00	Purpose and use	6
1.01	Manufacturer	6
1.02	Intended purpose	6
1.03	Area of use	6
1.04	Conditions of use	6
1.05	Important/Precautions	6
1.06	UHF RFID tag	7
1.07	Labels and Marking	7
1.08	Use	7
2.00	Maintenance	7
2.01	Cleaning	7
2.02	The owner's daily maintenance duty	7
2.03	Disposal of slings	8
3.00	Service and lifetime	8
3.01	Safety/service inspections	8
3.02	Lifetime	8
4.00	Technical specifications	9
5.00	EU-Declaration of conformity	9
6.00	Environmental policy statement	
	– V. Guldmann A/S	9
7.00	Warranty and service conditions	9
A.	Warranty	9
B.	Service or Repair	9
8.00	Placing the sling	34
9.00	Placing the sling	
	– Hip Support Sling with the Standing Sling	42
10.00	Product combinations	51

1.00 Purpose and use**1.01 Manufacturer**

V. Guldmann A/S
 Graham Bells Vej 21-23A
 DK-8200 Aarhus N
 Tel. + 45 8741 3100
 www.guldmann.com

1.02 Intended purpose

The sling is intended for lifting or supporting a person or body parts of a person.

1.03 Area of use

The sling is suited for use in hospitals, nursing homes, institutions, rehabilitation centers and in private homes.

The sling is used in combination with Guldmann GLS5 active lifter to assist the user from a sitting to a standing position.

1.04 Conditions of use

The sling and the active lifter are designed to lift or transfer a person who is capable of actively participating in the process. The user must also have good communication skills.

The use of the sling is subject to the following:

- The user can balance in a sitting position and has strength in the lower and upper body and can communicate
- The sling is used by trained staff or persons who have been instructed in the use of the sling in question.
- The correct size of sling is used.
- The maximum nominal load of 205 kg must not be exceeded.
- The helper pays attention to the well-being of the user when using the sling.
- The sling is used with the active lifter from Guldmann.

Important!

Plan the transfer. Avoid leaving the user in the sling unattended. Do not start to lift until it has been checked that the user cannot get trapped and that the sling does not catch on the bed, wheel-chair etc. The user's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes and wires that are attached to the user and/or equipment. Check that the hand control and hand control cable is free of the lifter, user and other objects before the lifting boom is activated up or down moved.

Guldmann shall not be liable for faults or accidents due to incorrect use of the lifting sling, or for reasons of inadequate attention on the part of the carer or user. If the sling is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.


1.05 Important/Precautions

- Read the instructions carefully before using the sling.
- The sling's maximum load must never be exceeded.
- The sling may only be used to lift a person.
- Before a sling is used, it must be examined according to point 2.02.
- Never use a sling that is too big for the user.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

1.06 UHF RFID tag

This product features a passive UHF RFID tag integrated in the product label. The RFID tag can be used for asset management & tracking purposes. The RFID tag is readable with equipment compliant to EPC global UHF Class I Gen ISO 18000-63.

1.07 Labels and Marking

CE marking

Medical Device Class I in accordance with EU MDR Regulation



Read the manual before use


UK Responsible Person
European Device Solutions Ltd. 15 Coanwood Drive, Whitley Bay, Tyne & Wear, NE25 9GB, United Kingdom. Email: info@europeandevicesolutions.co.uk
Tel: +44-754-559-5464

Example of serial number label Inspection label

Product label

Size label

1.08 Use

If there is any doubt about the selection or use of a lifting sling, please contact your supplier.

Caution!

Be careful when attaching the sling's lifting straps on the hooks of the lifting boom. Check that the straps have been correctly placed on the hooks. When pressing the up button on the hand control to lift the user, check again that all straps remain correctly placed in the lifting boom's hooks.

Applying the sling, see page 34.

2.00 Maintenance

2.01 Cleaning

Normal washing at the indicated temperature

Do not use bleaching agent

Tumble-drying at low temperature

Do not iron

2.02 The owner's daily maintenance duty

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Sling inspection checklist

Before using a Guldmann sling / accessory check the following:

Is the sling clean?

Follow facility specific infection control procedure.

Is the sling's label present, legible and complete?

Missing, illegible or incomplete sling label(s) could make identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling impossible.

Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?

Conclusion

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

2.03 Disposal of slings

Slings are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

3.00 Service and lifetime

3.01 Safety/service inspections

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection **must** be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Safe Operating Practices with Slings

Considerations for damaged or defective slings and taking them out of service:

Withdraw the sling from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

Sling inspection system

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

1) Initial

This level of inspection is done at the time that the sling is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the sling work load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage.

2) Frequent

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.

3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling.

Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

Sling inspection technique

The sling inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.


Consider the practice of documenting sling inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the unique sling identification number (important in differentiating similar slings), as well as the condition of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

Sample visual examples of synthetic sling damage *)


Chemical/caustic burns

Broken stitching


Crushed / Frayed webbing

Knots

Melting / Charring

*) sample visual images not intended to represent all types of potential damage

3.02 Lifetime

The expected lifetime of the sling is 5 years, however individually depending on usage pattern, washing etc. Before use the sling must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded if necessary.

Operation

The products operational environment:

- Operation temperatures between +10°C and +35°C (50°F and 95°F)
- A relative air humidity of between 30% and 70%

Beside temperature, the same environmental conditions apply for transportation and storage. The sling shall be stored on a flat, clean surface or hung on hooks using the lifting loops.

- Transport and storage temperatures between -10°C and +40°C (14°F and 104°F)

4.00 Technical specifications

Lifting capacity, SWL 205 kg (560 lbs)
Material Polyester^{x)}

^{x)} Fire retardant according to EN 1021

5.00 EU-Declaration of conformity

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

6.00 Environmental policy statement – V. Guldmann A/S

At Guldmann we will work actively to ensure that the negative impact that we can control is minimised.

Guldmann's Ambition is to ensure ongoing improvement of our environmental management system and its performance by:

- Working closely with our suppliers to ensure that we use materials and processes that are as sustainable as possible
- Continuously minimising the relative amount of waste and emissions and to ensure the highest possible degree of recycling
- Ensuring that our products do not have an unnecessary negative environmental impact in connection with use, recirculation and possibly destruction
- Complying with the applicable legislation
- Ensuring ongoing improvement of our environmental management system and associated environmental performance

All subsidiaries in the Guldmann group are covered by the above policy, and we expect that our Partners (suppliers and distributors) live up to this policy.

All Guldmann employees are obliged to immediately inform the management if they become aware of any violation of the environmental policy internally in the organisation or at our Partners.

This considers the economic and technological resources at our disposal and our general financial goals for the company and based on our fundamental values.

7.00 Warranty and service conditions

A. Warranty

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been

altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

Only for USA

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann's Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

B. Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number. If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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4. FR

Approchez le harnais Active Lifter.

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
5. GB/US

Place the user's feet onto the foot plate and adjust the knee pad - the top of the knee pad should be level with the top of the shin on a distance of two fingers below the knee joint.

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
6. GB/US

Optional – if the user has a tendency to have involuntary movements you can use the calf strap to ensure the user's feet do not come off the foot plate. If the calf strap is not in use fasten the strap to the back of the knee pad (avoiding loose straps).

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
7. GB/US

Raise the lifting boom to gain tension on the sling – **Ensure that the sling is still secure and attached before continuing with the lift.**

Communication with the user should be maintained throughout the lifting procedure.

Alternatively you have the option to attach the sling onto the lower hooks; this will facilitate an alternative lifting position and is more suitable for the shorter user.

When placing the user back into a seated position push the lifter up to the chair / commode chair / toilet and position the user so the user's calf is parallel to the chair / commode chair / toilet. Begin to lower and ask the user to push his/her pelvis backwards to ensure the correct seated position.

To remove the sling carry out the fitting instructions in reverse order, taking care to remove the sling from the user avoiding skin friction and movement.

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Stroppen er nu klar til brug.

1. GB/US

Attachment of the hip support straps

The support straps have either one or two loops at one end and a buckle at the other. Guide the end with the loop (s) through the hoop of the lifting arm and secure the support strap over the push handle as shown in the illustration.

Fig. 1 GLS5.1

On model GLS5.1 hip support straps with 2 loops should be used.

When mounting the hip support straps on the push handle, use the loop located 250 mm from the end as shown in the illustration.

Fig. 2 GLS5

On model GLS5 use the loop located at the end of the strap.

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4. GB/US

Optional – if the user has a tendency to
have involuntary movements you can use
the calf strap to ensure the user's feet do
not come off the foot plate. If the calf strap
is not in use fasten the strap to the back of
the knee pad (avoiding loose straps).

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

5. GB/US

Attach the hip support sling to the straps snapping the buckle together. Raise the lifting boom to gain tension on both slings – **Ensure that the sling is secure and attached before continuing with the lift.**

Adjust the hip support straps through the buckle and continue. The hip support sling should tension before the standing sling, – if it does not, then lower the user and readjust one or both slings.

As the lift progresses at approximately 60% into the standing procedure the hip support sling will loosen and once in standing the sling will be redundant allowing access to lower garments. Communication should be maintained throughout the lifting procedure.

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6. GB/US

When placing the user back into a seated position push the lifter up to the chair/commode chair/toilet and position the user so his/her calf is parallel to the chair/commode chair/toilet.

Begin to lower and at approximately 40% into the lowering procedure re-apply the hip support sling and continue to lower. The sling will ensure correct seated position.

To remove the sling carry out the fitting instructions in reverse order, taking care to peel the sling away from the user avoiding skin friction and movement.

Alternatively: The Active Micro Plus sling can be used with the active lifter – for instructions see separate user guide for Active Micro Plus Sling.

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Basic sling, Polyester	2700x1
Basic Low sling, Polyester	2710x1
Basic High sling, Polyester	2720x1
Basic Hammock, sling	2740x1
Basic sling, Net, fixed padding	2701x3
Basic Low sling, Net, fixed padding	2711x3
Basic High sling, Net, fixed padding	2721x3
Basic Shell, sling	2750x2
Basic Comfort High, polyester	2770x1
Basic Comfort High, net	2770x2

Custom Amputee seji	2900x1
Sit-On Comfort	2930x1
Sit-On Comfort High	2940x1
Sit-On II	2970x1
Sit-On High II	2980x1
Sit-On Comfort High, hygiene	2941x1
Sit-On	2950x1
Sit-On High	2960x1
Modified Sling	2949x

Repo. Sling, Bariatric	284656
Repo. Sling	28465
Repo. Sling, Short,	284653
Repo. Sling, Poly	284660
Repo. Sling, Grey net	284651
Repo. Sling, Grey net	284658
Repo. Sling, Grey net	284662
Repo. Sling, TENCEL	284657
Repo. Sling, Spacer	284659
Repo. Sling, Spacer, 6 loops	284669

Horizontal Sling, Standard	28463
Lifting sheet	2844851
Multi Support Sling, one size	28467
OR Sling, Poly	2848x1
Leg Sling Box of 10 pcs	28650
Pannus Support	28660
Turner	28700
Twin Turner	28751
Twin Turner, Bariatric	28760

Disposable High + Kids	2836x5
Disposable Twin Turner II, regular	287501
Disposable Twin Turner II, large	287511
Disposable Leg sling II	286501
Disposable High, bariatric	2836x2
Disposable Repositioning sling, 500 kg	284555
Disp. Horizontal sling, standard, 350 kg	284631
Disposable Multi support sling	284223
Disposable OR Sling	2848x5
Disposable Gait Trainer	2835x5
Disposable Comfort High	2770x5
Disposable Micro Plus	2815x5
Prone Positioning Sling (284221)	284225
Side Positioning Sling (Kit)	284226

Hanger						
Lifting hanger X-SMALL Item no. 556870	Lifting hanger SMALL Item no. 556880	Lifting hanger MEDIUM Item no. 556890	H-hanger Item no. 556950	Cross hanger 400 kg Item no. 561610	Cross hanger 500 kg Item no. 550800	Connecting bar + Cross hanger 500 kg Item no. 550544
x	x	x	x	x	x	
						x

x	x	x	x	x		
---	---	---	---	---	--	--

2-4, 4-6, 6-10, 10-16, XS	10-16, XS, S, M	S, M, L, XL, 2XL (6-10, 10-16)	S, M, L, XL - 5XL	M, L, XL - 5XL	M, L, XL - 5XL	
x	x	x	x	x	x	
x	x	x	x	x	x	x
x	x	x (x)	x	x	x	x
x	x	x (x)				
		x	x	x	x	x
	x	x	x	x	x	x
x	x	x				

x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x

x	x	x				
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x

		x	x	x	x	x
		x	x	x	x	x
	x	x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x

x	x	x		x	x	x
				x		
x	x	x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x

x	x	x	x	x		
		x		x	x	
		x		x	x	
		x		x	x	
		x	x	x	x	x
		x	x	x	x	x
x	x	x		x	x	
				x		
	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x	x	x
		x		x	x	x
		x		x	x	x

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Guldmann GmbH
Guldmann Sverige AB
Guldmann Srl
Guldmann Inc.