
DK .....	Positionerings sejl .....	2
GB/US .....	Positioning Sling .....	6
DE.....	Positionierungstuch .....	10
SE.....	Positionerings Sele .....	14
NO .....	Posisjoneringsseil .....	18
FR.....	Harnais de positionnement .....	22
IT .....	Imbragatura di posizionamento .....	26
NL.....	Positioning Sling .....	30

User manual – vers. 103.0

**DK... POSITIONERINGSSEJL**

Vers. 103.0

**Sidelejring**

1 x Positionerings sejl Varenr. : 284226

**Buglejring**

Pakke for bugleje Varenr. : 284225  
 2 x *Positionerings sejl*  
 2 x *Multisupport sejl*  
 1 x *Stabiliserings strop*

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<b>1.00</b>	<b>Purpose and use</b> . . . . .	<b>6</b>
1.01	Manufacturer . . . . .	6
1.02	Intended purpose . . . . .	6
1.03	Area of use . . . . .	6
1.04	Conditions of use . . . . .	6
1.05	Important/Precautions . . . . .	7
1.06	Labels and Marking . . . . .	7
1.07	Use . . . . .	7
<b>2.00</b>	<b>Maintenance</b> . . . . .	<b>7</b>
2.01	Cleaning . . . . .	7
2.02	The owner's daily maintenance duty . . . . .	7
2.03	Disposal of slings . . . . .	8
<b>3.00</b>	<b>Service and lifetime</b> . . . . .	<b>8</b>
3.01	Safety/service inspections . . . . .	8
3.02	Lifetime . . . . .	9
<b>4.00</b>	<b>Technical specifications</b> . . . . .	<b>9</b>
<b>5.00</b>	<b>EU-Declaration of conformity</b> . . . . .	<b>9</b>
<b>6.00</b>	<b>Environmental policy statement</b> – V. Guldmann A/S . . . . .	<b>9</b>
<b>7.00</b>	<b>Warranty and service conditions</b> . . . . .	<b>9</b>
A.	Warranty . . . . .	9
B.	Service or Repair . . . . .	9
<b>8.00</b>	<b>For instructions on how to use the positioning sling . to move a patient sideways.</b> . . . . .	<b>34</b>
<b>9.00</b>	<b>For instructions on how to the use the positioning . . sling to move a patient to a face-down position . . .</b> . . . . .	<b>44</b>
<b>10.00</b>	<b>Product combinations</b> . . . . .	<b>59</b>

**1.00 Purpose and use****1.01 Manufacturer**

V. Guldmann A/S  
 Graham Bells Vej 21-23A  
 DK-8200 Aarhus N  
 Tel. + 45 8741 3100  
 www.guldmann.com

**1.02 Intended purpose**

The sling is intended for lifting or supporting a person or body parts of a person

**1.03 Area of use**

The sling is suited for use in hospitals.

**1.04 Conditions of use**

The sling is designed for use with a ceiling hoist systems. The positioning sling is used to turn and position patients at the OR table before and after surgical procedures.

The use of the sling is subject to the following:

- The sling is used by trained staff or persons who have been instructed in the use of the sling in question.
- The safe working load of 255 kg must never be exceeded.
- The sling is used to turn a person in a lying position.
- The helper pays attention to the well-being of the user when using the sling.
- The sling is used with the Guldmann lifting hanger.

**Side positioning**

During side positioning procedures, patients are turned using a positioning sling and moved from lying on their back to lying on their side.

**Prone positioning**

During prone positioning procedures use the package that contains two positioning slings, two Multi support slings and one security strap. The patient is turned from a supine to a prone position by using one positioning sling and then back to a supine position using the other positioning sling following the operation.

Once the patient has been turned over and is in the prone position, it is possible to lift parts of the body using the accompanying Multi support sling, so that positioning pillows can be placed underneath the patient prior to the operation.

**Important!**

Plan the move. Never leaving the user in the lifting sling unattended. Do not start to lift until it has been checked that the user cannot get trapped and that the sling does not catch on the bed, wheelchair or other obstacles. The user's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes and wires that are attached to the user and/or equipment. Check that the hand control and hand control cable is free of hanger, patient and other object before the hoist is activated up or down moved.

Guldmann shall not be liable for faults or accidents due to incorrect use of the lifting sling, or for reasons of inadequate attention on the part of the carer or user. If the sling is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

## 1.05 Important/Precautions

- Read the instructions carefully before using the sling.
- The sling's maximum load must never be exceeded.
- The sling may only be used to roll or support a person.
- Before a sling is used, it must be examined according to point 2.02.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

## 1.06 Labels and Marking

CE marking

Medical Device Class I in accordance with EU MDR Regulation


Read the manual before use


**UK Responsible Person**  
**European Device Solutions Ltd.** 15 Coanwood Drive,  
 Whitley Bay, Tyne & Wear, NE25 9GB, United Kingdom.  
 Email: info@europeandevicesolutions.co.uk  
 Tel: +44-754-559-5464

### Product label

### Example of LOT number label

## 1.07 Use


**If there is any doubt about the selection or use of a lifting sling, please contact your supplier.**

### Lifting hanger, 4 attachment points

#### Caution!

**Be careful when attaching the lifting sling's straps on the hooks. Check that the straps have been correctly placed in the lifting hanger's hooks. When pressing the up button on**


**the hand control to lift the user, check again that all straps remain correctly placed in the lifting hanger's hooks (Fig. 1).**

### Lifting hanger

#### Caution!

**Be careful when attaching the lifting sling on the hooks. Check that the straps have been pulled completely through the rubber safety catch (A) and into place in the lifting hanger's hooks. When pressing the up button to lift the user, check again that all the straps remain correctly placed in the lifting hanger's hooks (fig. 1a and fig. 1b).**

*For instructions on how to use the positioning sling for side positioning procedures, see page 34.*

*For instructions on how to use the positioning sling for prone positioning procedures, see page 44.*

## 2.00 Maintenance

### 2.01 Cleaning

Normal washing at the indicated temperature  
 The sling can be washed a maximum of 5 times

Do not use bleaching agent

Tumble-drying at low temperature

Do not iron

### 2.02 The owner's daily maintenance duty

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

#### Sling inspection checklist

Before using a Guldmann sling / accessory check the following:

#### Is the sling clean?

Follow facility specific infection control procedure.

#### Is the sling's label present, legible and complete?

Missing, illegible or incomplete sling label(s) could make identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling impossible.

#### Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

**Is the fabric intact?**

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

**Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?****Conclusion**

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

**2.03 Disposal of slings**

Slings are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

**3.00 Service and lifetime****3.01 Safety/service inspections**

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection must be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

**Safe Operating Practices with Slings**

Considerations for damaged or defective slings and taking them out of service:

Withdraw the sling from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

**Sling inspection system**

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

**1) Initial**

This level of inspection is done at the time that the sling is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the sling work load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage.

**2) Frequent**

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.


**3) Periodic**

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

**Sling inspection technique**

The sling inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the unique sling identification number (important in differentiating similar slings), as well as the condition of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

**Sample visual examples of synthetic sling damage \*)****Chemical/caustic burns****Broken stitching****Crushed / Frayed webbing**

Knots

Melting / Charring

*x) sample visual images not intended to represent all types of potential damage*

**3.02 Lifetime**

The life of the sling is individual and depends on how it is used, washed etc. Before use the sling must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded.

The sling can be washed a maximum of 5 times. Remember to mark off the label following each wash. The sling must be disposed of after 5 washes.

**4.00 Technical specifications**

Lifting capacity, SWL . . . . . 255 kg  
 Material, white part: . . . . . Polyester  
 Sliding material . . . . . Nylon

**5.00 EU-Declaration of conformity**

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

**6.00 Environmental policy statement  
 – V. Guldmann A/S**

At Guldmann we will work actively to ensure that the negative impact that we can control is minimised.

Guldmann’s Ambition is to ensure ongoing improvement of our environmental management system and its performance by:

- Working closely with our suppliers to ensure that we use materials and processes that are as sustainable as possible
- Continuously minimising the relative amount of waste and emissions and to ensure the highest possible degree of recycling
- Ensuring that our products do not have an unnecessary negative environmental impact in connection with use, recirculation and possibly destruction
- Complying with the applicable legislation
- Ensuring ongoing improvement of our environmental management system and associated environmental performance

All subsidiaries in the Guldmann group are covered by the above policy, and we expect that our Partners (suppliers and distributors) live up to this policy.

All Guldmann employees are obliged to immediately inform the management if they become aware of any violation of the environmental policy internally in the organisation or at our Partners.

This considers the economic and technological resources at our disposal and our general financial goals for the company and based on our fundamental values.

**7.00 Warranty and service conditions**

**A. Warranty**

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the “Warranty Period”). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

**Only for USA**

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann’s Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

**B. Service or Repair**

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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## 1. GB/US

**Preparation for side positioning**

If a positioning bolster will be used during procedure, place the Multi Support sling **under** the OR Turning sling **before** it is placed on the OR table.

The Multi support sling must be placed at the chest and below the axilla. (see photo 9 on page 35 for proper positioning).

**Important!**

The center mark must always be face-up and positioned in the middle. White and soft side face upwards. The positioning sling's side border with the long folded white and blue edge should be on the side the patient is turning **toward** and will lay on.

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## 2. GB/US

Attach the positioning sling to the OR table by fixing the sling's 4 elastic bands under the 4 corners of the OR table.

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## 3. GB/US

**Important!**

Always utilize the necessary safety steps to accompany use of the sling, such as correct number of caregivers and use of OR table side stops.

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## 4. GB/US

**Side positioning – turning patient from a supine to a side-lying position**

With the patient lying supine (face up) on the positioning sling.

Prepare the patient's body for the side-lying position, taking care to position the dependent upper and lower extremity for side-lying.

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## 5. GB/US

The positioning sling's lifting straps are now placed onto the lifting hanger. Ensure that the patient's "anaesthesia arm" is on the same side as the lifting straps. By doing so, you can ensure that this arm is always on top during the turning procedure, so that tubes and cables do not become tangled or caught. Loosen the 4 elastic bands from the corners.

**Important!**

Remember to check the elastic bands to ensure they have been loosened from all 4 corners.

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## 6. GB/US

Push on the patient's hip and shoulders while simultaneously activating the lift, so that the lifting hanger and the attached positioning sling are moved upwards. The turning procedure has now begun.

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## 7. GB/US

Once in sidelying position, the top leg can be moved onto a bolster or other.

Should it become necessary to lift the patient's upper body, a Multi support sling must be placed under the positioning sling on the OR table beforehand in order to lift the patient's upper body (see photo 9 on page 35).

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#### 8. GB/US

When the patient has been turned over and is lying on his/her side, manually remove the sling by pulling it out piece by piece, material towards material. Avoid creating friction against the patient's skin.

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## 9. GB/US

If necessary: Lift the patient's upper body with the Multi support sling, and place a bolster, pad or other device as needed. Lower the patient back onto the Multi support sling during the operation. Remove the straps from the lifting hanger.

Ensure that the patient is lying correctly and comfortably, and without any entanglements.

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## 2. GB/US

Position and fasten the bottom part of the two-piece positioning sling to the OR table by fixing the positioning sling's 4 elastic bands around the 4 corners of the OR table.

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## 3. GB/US

Should it become necessary to lift the patient's upper body when placing bolsters, place the two Multi support slings over the bottom positioning sling at the patient's upper body so that they support the pelvis and abdomen, as well as at the patient's chest below the axilla (see photo 9 on page 46 for proper positioning).

The patient's arms must be outside of the sling.

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## 4. GB/US

Attach the top half of the two-piece positioning sling in the same way as the bottom half by fixing the 4 elastic bands around the OR tables 4 corners. Note however that the two lifting straps must now be placed on the opposite side.

**Important!**

Always utilize the necessary safety steps to accompany use of the sling, such as correct number of caregivers and use of OR table side stops. There must be mounted sidestops on each side of the OR table.

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## 6. GB/US

When the patient is ready for turning, the upper positioning sling's lifting straps are now placed onto the ceiling hoist's hanger. Ensure that the patient's "anaesthesia arm" is on the same side as the lifting straps. By doing so, you can ensure that this arm is always on top during the turning procedure, so that tubes and cables do not become tangled or caught. Loosen the 4 elastic bands from the corners.

**Important!**

Remember to ensure that only the elastic bands belonging to the top sling are loosened. These 4 bands are all of the same colour, though their colour is different from the bands on the bottom sling (the top sling's elastic bands are either white or black, depending on which side the turning procedure is being done on).

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## 7. GB/US

**Important!**

Do not push the patient before the "side stoppers" have been fitted onto the OR table. Push the patient's hip and shoulders while simultaneously activating the lift, so that the hanger and the mounted positioning sling are moved upwards. The turning procedure has now begun.

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#### 8. GB/US

When the patient has been turned over and is lying prone position, manually remove the sling by pulling it out piece by piece, material towards material. Avoid creating friction against the patient's skin.

#### 8. DE

When the patient has been turned over and is lying prone position, manually remove the sling by pulling it out piece by piece, material towards material. Avoid creating friction against the patient's skin.

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## 9. GB/US

If necessary, the patient's arms can be placed over his/her head. Adjust the two Multi support slings so that they are optimally positioned. Place the straps onto the lifting hanger. Be aware that the hanger is placed parallel with the patient's body.

Lift the patient's upper body and place the bolster / bolsters or other device. Lower and place the patient as required. The Multi support sling stays under the patient during the procedure. Remove the straps from the lifting hanger.

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## 11. GB/US

### Prone positioning – turning the patient from a prone to a supine position

Remove the security strap before moving the patient back to the supine position. Loosen the 4 elastic bands from the corners. Reattach the Multi support slings to the lifting hanger, lift the patient, and remove the bolsters. Lower the patient onto the bottom half of the two positioning slings.

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## 13. GB/US

Prepare the patient for turning supine. The patient can be turned by pushing again on his/her hip and shoulders while simultaneously activating the lift.

**Important! "Side stoppers" must be fitted onto the OR table.**

Manually pull out the positioning slings and the two Multi support slings. Take care to avoid causing friction against the patient's tissue. Ensure that the patient is lying correctly and comfortably, and without any entanglements.

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Basic sling, Polyester	2700x1
Basic Low sling, Polyester	2710x1
Basic High sling, Polyester	2720x1
Basic Hammock, sling	2740x1
Basic sling, Net, fixed padding	2701x3
Basic Low sling, Net, fixed padding	2711x3
Basic High sling, Net, fixed padding	2721x3
Basic Shell, sling	2750x2
Basic Comfort High, polyester	2770x1
Basic Comfort High, net	2770x2

Custom Amputee seji	2900x1
Sit-On Comfort	2930x1
Sit-On Comfort High	2940x1
Sit-On II	2970x1
Sit-On High II	2980x1
Sit-On Comfort High, hygiene	2941x1
Sit-On	2950x1
Sit-On High	2960x1
Modified Sling	2949x

Repo. Sling, Bariatric	284656
Repo. Sling	28465
Repo. Sling, Short,	284653
Repo. Sling, Poly	284660
Repo. Sling, Grey net	284651
Repo. Sling, Grey net	284658
Repo. Sling, Grey net	284662
Repo. Sling, TENCEL	284657
Repo. Sling, Spacer	284659
Repo. Sling, Spacer, 6 loops	284669

Horizontal Sling, Standard	28463
Lifting sheet	2844851
Multi Support Sling, one size	28467
OR Sling, Poly	2848x1
Leg Sling Box of 10 pcs	28650
Pannus Support	28660
Turner	28700
Twin Turner	28751
Twin Turner, Bariatric	28760

Disposable High + Kids	2836x5
Disposable Twin Turner II, regular	287501
Disposable Twin Turner II, large	287511
Disposable Leg sling II	286501
Disposable High, bariatric	2836x2
Disposable Repositioning sling, 500 kg	284555
Disp. Horizontal sling, standard, 350 kg	284631
Disposable Multi support sling	284223
Disposable OR Sling	2848x5
Disposable Gait Trainer	2835x5
Disposable Comfort High	2770x5
Disposable Micro Plus	2815x5
Prone Positioning Sling (284221)	284225
Side Positioning Sling (Kit)	284226

Hanger						
Lifting hanger X-SMALL Item no. 556870	Lifting hanger SMALL Item no. 556880	Lifting hanger MEDIUM Item no. 556890	H-hanger Item no. 556950	Cross hanger 400 kg Item no. 561610	Cross hanger 500 kg Item no. 550800	Connecting bar + Cross hanger 500 kg Item no. 550544
						
x	x	x	x	x	x	
						x

x	x	x	x	x		
---	---	---	---	---	--	--

2-4, 4-6, 6-10, 10-16, XS	10-16, XS, S, M	S, M, L, XL, 2XL (6-10, 10-16)	S, M, L, XL - 5XL	M, L, XL - 5XL	M, L, XL - 5XL	
x	x	x	x	x	x	
x	x	x	x	x	x	x
x	x	x (x)	x	x	x	x
x	x	x (x)				
		x	x	x	x	x
	x	x	x	x	x	x
x	x	x	x	x		

x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x		

x	x	x				
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x

		x	x	x	x	x
		x	x	x	x	x
	x	x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x

x	x	x		x	x	x
				x		
x	x	x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x

x	x	x	x	x		
		x		x	x	
		x		x	x	
		x		x	x	
		x	x	x	x	x
		x	x	x	x	x
x	x	x		x	x	
				x		
	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x	x	x
		x		x	x	x
		x		x	x	x

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**Guldmann GmbH**
**Guldmann Sverige AB**
**Guldmann Srl**
**Guldmann Inc.**