DK	Engangssejl Gait Trainer	3
GB/US	Disposable Gait Trainer	7
DE	Einweg Gait Trainer	12
SE	Engångsele Gait Trainer	16
NO	Engangs gåseil	20
FR	Gait trainer jetable	24
IT	Gait trainer monouso	28
NL	Disposable Gait Trainer	32

User manual – vers. 103.0

Guldmann™

Vers. 103.0

Varenr.
 283545-1 S
 283555-1 M
 283565-1 L
 283575-1 XL
 283585-1 2XL

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Item nos:

283545-1 S
 283555-1 M
 283565-1 L
 283575-1 XL
 283585-1 2XL

1.00	Purpose and use	6
1.01	Manufacturer	6
1.02	Intended purpose	6
1.03	Area of use	6
1.04	Conditions of use	6
1.05	Important/Precautions	6
1.06	Labels and Marking	7
1.07	Use	7
2.00	Maintenance	7
2.01	Cleaning	7
2.02	The owner's daily maintenance duty	7
2.03	Disposal of slings	8
3.00	Service and lifetime	8
3.01	Safety/service inspections	8
3.02	Lifetime	9
4.00	Technical specifications	9
5.00	EU-Declaration of conformity	9
6.00	Environmental policy statement – V. Guldmann A/S	9
7.00	Warranty and service conditions	9
A.	Warranty	9
B.	Service or Repair	9
8.00	Placing the sling	34
9.00	Product combinations	45

1.00 Purpose and use

1.01 Manufacturer

V. Guldmann A/S
 Graham Bells Vej 21-23A
 DK-8200 Aarhus N
 Tel. + 45 8741 3100
 www.guldmann.com

1.02 Intended purpose

The sling is intended for lifting or supporting a person or body parts of a person.

1.03 Area of use

The sling is suited for use in hospital, a nursing home, institutions and in rehabilitation centers.

1.04 Conditions of use

As the sling is a disposable sling, it is suitable as a personal sling and in situations where a high degree of hygiene is required and for infection control program. The user's name can be written on the sling with the included pen. If necessary the sling can be discarded every time it has been used – or when the user is discharged.

The sling is designed to be used with ceiling hoist systems and is ideal for use in connection with standing and walking training.

The sling is designed to lift/ support users that suffer from reduced balance but are able to stand.

The use of the sling is subject to the following:

- The sling is used by trained staff or persons who have been instructed in the use of the sling in question.
- The correct size of sling is used.
- The maximum nominal load, 205 kg (450 lbs) must not be exceeded.
- The sling is used for lifting/ supporting in connection with standing and walking training.
- The helper pays attention to the well-being of the user when using the sling.
- The sling is used with the Guldmann lifting hanger.

Important!

Plan the move. Never leaving the user in the lifting sling unattended. Do not start to lift until it has been checked that the user cannot get trapped and that the sling does not catch on the bed, wheelchair or other obstacles. The user's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes and wires that are attached to the user and/or equipment. Check that the hand control and hand control cable is free of hanger, patient and other object before the hoist is activated up or down moved.

Guldmann shall not be liable for faults or accidents due to incorrect use of the lifting sling, or for reasons of inadequate attention on the part of the carer or user. If the sling is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

1.05 Important/Precautions

- Please read the instructions carefully before using the sling.
- The slings maximum load must never be exceeded.
- The sling may only be used to lift a person.

- Disposable slings must not be used for bath and in swimming pools.
- Before a sling is used, it must be examined according to point 2.02.
- Never use a sling that is too big for the user.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

1.06 Labels and Marking

CE marking

Medical Device Class I in accordance with EU MDR Regulation

Read the manual before use


Single Patient Multiple use

UK CA UK Responsible Person
European Device Solutions Ltd. 15 Coanwood Drive,
 Whitley Bay, Tyne & Wear, NE25 9GB, United Kingdom.
 Email: info@europeandevicesolutions.co.uk
 Tel: +44-754-559-5464

Product label

LOT number label


1.07 Use

If there is any doubt about the selection or use of a lifting sling, please contact your supplier.

Lifting hanger, 4 attachment points

Caution!


Be careful when attaching the lifting sling's straps on the hooks. Check that the straps have been correctly placed in the lifting hanger's hooks. When pressing the up button on the hand control to lift the user, check again that all straps remain correctly placed in the lifting hanger's hooks (Fig. 1).

Lifting hanger

Caution!

Be careful when attaching the lifting sling on the hooks. Check that the straps have been pulled completely through the rubber safety catch (A) and into place in the lifting hanger's hooks. When pressing the up button to lift the user, check again that all the straps remain correctly placed in the lifting hanger's hooks (fig. 1a and fig. 1b).

Placing the sling, look at page 34

2.00 Maintenance

2.01 Cleaning

- Do not wash
- Do not use bleaching agent
- Do not tumble dry
- Do not iron

The sling is a disposable sling and if necessary it can be discarded every time it has been used – or when the user is discharged.

The label "Do not wash" will be changed to "Do not use" if the sling has been washed.

2.02 The owner's daily maintenance duty

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Sling inspection checklist

Before using a Guldmann sling / accessory check the following:

Is the sling clean?

Follow facility specific infection control procedure.

Is the sling's label present, legible and complete?

Missing, illegible or incomplete sling label(s) could compromise identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling.

Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?**Conclusion**

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

2.03 Disposal of slings

Slings are disposed of by incineration.

By proper incineration polyester will be degraded to carbon dioxide and water.

3.00 Service and lifetime**3.01 Safety/service inspections**

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection **must** be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Safe Operating Practices with Slings

Considerations for damaged or defective slings and taking them out of service. Withdraw the sling from service if any of the following conditions exist:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

Sling inspection system

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

1) Initial

This level of inspection is done at the time that the sling is received into your facility. The inspector should insure that no damage has occurred during transit, and also verify that the sling work load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage

2) Frequent

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.


3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.


Sling inspection technique

The sling inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the unique sling identification number (important in differentiating similar slings), as well as the condition of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

Sample visual examples of synthetic sling damage *)**Chemical/caustic burns****Broken stitching**

Crushed / Frayed webbing

Knots

Melting / Charring

x) sample visual images not intended to represent all types of potential damage

equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

3.02 Lifetime

The sling is designed for short-term use by one particular user and is to be discarded when soiled or not longer needed.

4.00 Technical specifications

Lifting capacity, SWL 205 kg (450 lbs)
Material.....Polyester

Only for USA

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann's Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

5.00 EU-Declaration of conformity

CE The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

**6.00 Environmental policy statement
- V. Guldmann A/S**

Guldmann is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

It is Guldmann's goal to:

- Comply with the current environmental legislation (e.g. WEEE and REACH directives)
- Ensure that we, at the widest possible range, use RoHS compliant materials and components
- Ensure that our products do not have an unnecessary negative impact on the environment regarding use, recirculation or disposal
- Ensure that our products contribute to a positive working environment in the places they are utilised

Inspections are made annually by the Department for Nature and Environment from the Municipality of Aarhus using the Danish Environmental Protection Act, section 42 as a reference.

7.00 Warranty and service conditions

A. Warranty

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or

B. Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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9. GB/US

The adjustable safety strap is now guided through the chest strap on the opposite side in the height that fits the patient.

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
Basic sling, Polyester	2700x1
Basic Low sling, Polyester	2710x1
Basic High sling, Polyester	2720x1
Basic Hammock, sling	2740x1
Basic sling, Net, fixed padding	2701x3
Basic Low sling, Net, fixed padding	2711x3
Basic High sling, Net, fixed padding	2721x3
Basic Shell, sling	2750x2
Basic Comfort High, polyester	2770x1
Basic Comfort High, net	2770x2

x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x		

Custom Amputee seji	2900x1
Sit-On Comfort	2930x1
Sit-On Comfort High	2940x1
Sit-On II	2970x1
Sit-On High II	2980x1
Sit-On Comfort High, hygiene	2941x1
Sit-On	2950x1
Sit-On High	2960x1
Modified Sling	2949x

x	x	x				
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x
x	x	x	x	x	x	x

Repo. Sling, Bariatric	284656
Repo. Sling	28465
Repo. Sling, Short,	284653
Repo. Sling, Poly	284660
Repo. Sling, Grey net	284651
Repo. Sling, Grey net	284658
Repo. Sling, Grey net	284662
Repo. Sling, TENCEL	284657
Repo. Sling, Spacer	284659
Repo. Sling, Spacer, 6 loops	284669

		x	x	x	x	x
		x	x	x	x	x
	x	x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x
		x	x	x	x	x

Horizontal Sling, Standard	28463
Lifting sheet	2844851
Multi Support Sling, one size	28467
OR Sling, Poly	2848x1
Leg Sling Box of 10 pcs	28650
Pannus Support	28660
Turner	28700
Twin Turner	28751
Twin Turner, Bariatric	28760

x	x	x		x	x	x
x	x	x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x
		x		x	x	x

Disposable High + Kids	2836x5
Disposable Twin Turner II, regular	287501
Disposable Twin Turner II, large	287511
Disposable Leg sling II	286501
Disposable High, bariatric	2836x2
Disposable Repositioning sling, 500 kg	284555
Disp. Horizontal sling, standard, 350 kg	284631
Disposable Multi support sling	284223
Disposable OR Sling	2848x5
Disposable Gait Trainer	2835x5
Disposable Comfort High	2770x5
Disposable Micro Plus	2815x5
Prone Positioning Sling (284221)	284225
Side Positioning Sling (Kit)	284226

x	x	x	x	x		
		x		x	x	
		x		x	x	
		x		x	x	
		x	x	x	x	x
		x	x	x	x	x
x	x	x		x	x	
				x		
	x	x	x	x	x	x
x	x	x	x	x		
x	x	x	x	x	x	x
		x		x	x	x
		x		x	x	x

<h1>Guldmann™</h1>	
Product combinations	
Lifting module / Mobile lifter	
GH1, GH1 F, GH1 Q, GHZ, GH3, GH3+ lifting module	
GH3 Twin lifting module	
GL5.2 Mobile lifter	

Stretcher			
Horizontal lifting support, stepless weight adjustment Item no. 28456	Horizontal lifting support, stepped weight adjustment Item no. 28466	Horizontal lifter, Foldable, 1 strap version Item no. 28445	Horizontal Lifter, Foldable, 2 strap version Item no. 28444
			
x	x	x	
	x		x
x	x	x	

Sling	Item no.
Active Micro Plus	2810x1
Active Micro, Poly	2840x1
Active Trainer	2830x1
Active Vest Kids	2831x1
Gait Trainer, Bariatric	283100
Gait Trainer	2832x1
Vest for Stand Shell	2835x1

Basic sling, Polyester	2700x1
Basic Low sling, Polyester	2710x1
Basic High sling, Polyester	2720x1
Basic Hammock, sling	2740x1
Basic sling, Net, fixed padding	2701x3
Basic Low sling, Net, fixed padding	2711x3
Basic High sling, Net, fixed padding	2721x3
Basic Shell, sling	2750x2
Basic Comfort High, polyester	2770x1
Basic Comfort High, net	2770x2

Custom Amputee sejl	2900x1
Sit-On Comfort	2930x1
Sit-On Comfort High	2940x1
Sit-On II	2970x1
Sit-On High II	2980x1
Sit-On Comfort High, hygiene	2941x1
Sit-On	2950x1
Sit-On High	2960x1
Modified Sling	2949x

Repo. Sling, Bariatric	284656
Repo. Sling	28465
Repo. Sling, Short	284653
Repo. Sling, Poly	284660
Repo. Sling, Grey net	284651
Repo. Sling, Grey net	284658
Repo. Sling, Grey net	284662
Repo. Sling, TENCEL	284657
Repo. Sling, Spacer	284659
Repo. Sling, Spacer, 6 loops	284669

Horizontal Sling, Standard	28463
Lifting sheet	2844851
Multi Support Sling, one size	28467
OR Sling, Poly	2848x1
Leg Sling Box of 10 pcs	28650
Pannus Support	28660
Turner	28700
Twin Turner	28751
Twin Turner, Bariatric	28760

x	x		
		x	x

Disposable High + Kids	2836x5
Disposable Twin Turner II, regular	287501
Disposable Twin Turner II, large	287511
Disposable Leg sling II	286501
Disposable High, bariatric	2836x2
Disposable Repositioning sling, 500 kg	284555
Disp. Horizontal sling, standard, 350 kg	284631
Disposable Multi support sling	284223
Disposable OR Sling	2848x5
Disposable Gait Trainer	2835x5
Disposable Comfort High	2770x5
Disposable Micro Plus	2815x5
Prone Positioning Sling (284221)	284225
Side Positioning Sling (Kit)	284226

x	x		

| Time to care |

V. Guldmann A/S

Head Office:

Tel. +45 8741 3100

info@guldmann.com

www.guldmann.com

Guldmann GmbH

Tel. +49 611 974 530

info@guldmann.de

www.guldmann.de

Guldmann Sverige AB

Tel. +46 0322 55290

info@guldmann.se

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Guldmann Norge

Tel. +47 906 00 130

smy@guldmann.com

www.guldmann.com

Guldmann Sarl

Tél. +33 145 54 78 36

france@guldmann.com

www.guldmann.fr

Guldmann Srl

Tel. +39 0521 660132

Cell. +39 340 9087107

italia@guldmann.com

www.guldmann.it

Guldmann BV

Tel. +31 053 428 30 90

info@ergocare.nl

www.guldmann.nl

Guldmann Inc.

Tel. 800 664 8834

Tel. 813 880 0619

info@guldmann.net

www.guldmann.net

Guldmann™