DK	Engangs Twin Turner II / Engangs Twin Turner Bariatric	2
GB/US	Disposable Twin Turner II / Disposable Twin Turner Bariatric	6
DE	Einweg Twin Turner II / Einweg Twin Turner Bariatric	10
SE	Engångssele Twin Turner II / Engångs Twin Turner Bariatric	14
NO	Engangs Twin Turner II / Engangs Twin Turner Bariatric	18
FR	Harnais Twin Turner II jetable / Harnais Twin Turner II jetable Bariatric	22
IT	Twin Turner II Monouso / Twin Turner monouso	26
NL	Disposable Twin Turner II / Disposable Twin Turner Bariatric	30
ES	Disposable Twin Turner II / Twin Turner Bariatric desechable	34

User manual – vers. 104.0

DK... ENGANGS TWIN TURNER II / ENGANGS TWIN TURNER BARIATRIC

Vers. 104.0

Varenr.
287501-1 Standard
287511-1 Large
287521-1 Bariatric

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287511-1 Large

287521-1 Bariatric

1.00	Purpose and use	6
1.01	Manufacturer	6
1.02	Intended purpose	6
1.03	Area of use	6
1.04	Conditions of use	6
1.05	Important/Precautions	6
1.06	Labels and Marking	7
1.07	Use	7
2.00	Maintenance	7
2.01	Cleaning	7
2.02	The owner's daily maintenance duty	7
2.03	Disposal of slings	8
3.00	Service and lifetime	8
3.01	Safety/service inspections	8
3.02	Lifetime	9
4.00	Technical specifications	9
5.00	EU-Declaration of conformity	9
6.00	Environmental policy statement – V. Guldmann A/S	9
7.00	Warranty and service conditions	9
A.	Warranty	9
B.	Service or Repair	9
8.00	Placing the sling	38
9.00	Product combinations	44

1.00 Purpose and use

1.01 Manufacturer

V. Guldmann A/S
Graham Bells Vej 21-23A
DK-8200 Aarhus N
Tel. + 45 8741 3100
www.guldmann.com

1.02 Intended purpose

The sling is intended for lifting or supporting a person or body parts of a person.

1.03 Area of use

The sling is suited for use in hospital, nursing homes, institutions and in rehabilitation centers.

1.04 Conditions of use

As the sling is a disposable sling, it is suitable as a personal sling and in situations where a high degree of hygiene is required and for infection control program. The user's name can be written on the sling with the included pen. If necessary the sling can be discarded every time it has been used – or when the user is discharged.

The sling is suited for turning and supporting people with disabilities in hospitals, nursing homes, rehabilitation centres and at home.

Disposable Twin Turner can be used with ceiling hoist system and it makes it possible to turn person on both sides without removing the sling.

The use of the sling is subject to the following:

- The sling is used by trained staff or persons who have been instructed in the use of the sling in question.
- The correct size of sling is used.
- The maximum nominal load, 205 kg / 500 kg (450 lbs / 1100 lbs) must not be exceeded.
- Disposable Twin Turner is used for turning a person on both sides in lying position.
- The helper pays attention to the well-being of the user when using the sling.
- The sling is designed for use with Guldmann lifting hangers, see the applicable combinations at the back of the manual. If the sling is used in combination with products that are not developed by Guldmann, then a risk assessment must be carried out by qualified personnel.

Important:

Guldmann accepts no responsibility for errors or accidents that may occur as a result of:

- incorrect use of the lifting sling
- combining other suppliers' products with Guldmann products.
- lack of attention from the helper or the user.

If the sling is used in combination with products that are not manufactured by Guldmann, then the same method of attaching the sling to the lifting hanger must be used and a risk assessment must be carried out by qualified personnel.

Plan the move. Never leave the user in the lifting sling unattended. Do not start to lift until it has been checked that the user cannot get trapped and that the sling does not catch on the bed, wheelchair or other obstacles. The user's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes and wires that are attached to the user and/or equipment. Check that the hand control and hand control cable is free of hanger, patient and other objects before the hoist is activated up or down moved.

1.05 Important/Precautions

- Read the instructions carefully before using the sling.
- The slings maximum load must never be exceeded.
- The sling may only be used to lift and support a person.
- Disposable slings must not be used for bath and in swimming pools.
- Before a sling is used, it must be examined according to point 2.02.
- Possible repairs must only be made by the manufacturer.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

1.06 Labels and Marking

CE marking

Medical Device Class I in accordance with EU MDR Regulation

Read the manual before use

Single Patient Multiple use


UK Responsible Person
European Device Solutions Ltd. 15 Coanwood Drive,
Whitley Bay, Tyne & Wear, NE25 9GB, United Kingdom.
Email: info@europeandevicesolutions.co.uk
Tel: +44-754-559-5464

Example of product label

LOT number label


1.07 Use

If there is any doubt about the selection or use of a lifting sling, please contact your supplier.

Lifting hanger, 4 attachment points

Caution!


Be careful when attaching the lifting sling's straps on the hooks. Check that the straps have been correctly placed in the lifting hanger's hooks. When pressing the up button on the hand control to lift the user, check again that all straps remain correctly placed in the lifting hanger's hooks (Fig. 1).

Lifting hanger

Caution!

Be careful when attaching the lifting sling on the hooks. Check that the straps have been pulled completely through the rubber safety catch (A) and into place in the lifting hanger's hooks. When pressing the up button to lift the user, check again that all the straps remain correctly placed in the lifting hanger's hooks (fig. 1a and fig. 1b).

Placing the sling, look at page 38.

2.00 Maintenance

2.01 Cleaning

- Do not wash
- Do not use bleaching agent
- Do not tumble dry
- Do not iron

The sling is a disposable sling and if necessary it can be discarded every time it has been used – or when the user is discharged.

The label "Do not wash" will be changed to "Do not use" if the sling has been washed.

2.02 The owner's daily maintenance duty

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Sling inspection checklist

Before using a Guldmann sling / accessory check the following:

Is the sling clean?

Follow facility specific infection control procedure.

Is the sling's label present, legible and complete?

Missing, illegible or incomplete sling label(s) could compromise identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling.

Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?

Conclusion

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

2.03 Disposal of slings

Slings are disposed of by incineration.

By proper incineration polypropylene and polyester will be degraded to carbon dioxide and water.

3.00 Service and lifetime**3.01 Safety/service inspections**

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection **must** be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Safe Operating Practices with Slings

Considerations for damaged or defective slings and taking them out of service. Withdraw the sling from service if any of the following conditions exist:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

Sling inspection system

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

1) Initial

This level of inspection is done at the time that the sling is received into your facility. The inspector should insure that no damage has occurred during transit, and also verify that the sling work

load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage

2) Frequent

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.


3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.


Sling inspection technique

The sling inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling inspections through written inspection records. The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the unique sling identification number (important in differentiating similar slings), as well as the condition of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

Sample visual examples of synthetic sling damage *)*Chemical/caustic burns**Broken stitching**Crushed / Frayed webbing**Knots*

x) sample visual images not intended to represent all types of potential damage

3.02 Lifetime

The sling is designed for short-term use by one particular user and is to be discarded when soiled or not longer needed.

4.00 Technical specifications

Lifting capacity, SWL 205 kg / 500 kg (450 lbs / 1100 lbs)
Material. Polypropylene/polyester

5.00 EU-Declaration of conformity

CE The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

6.00 Environmental policy statement
– V. Guldmann A/S

At Guldmann we will work actively to ensure that the negative impact that we can control is minimised.

Guldmann’s Ambition is to ensure ongoing improvement of our environmental management system and its performance by:

- Working closely with our suppliers to ensure that we use materials and processes that are as sustainable as possible
- Continuously minimising the relative amount of waste and emissions and to ensure the highest possible degree of recycling
- Ensuring that our products do not have an unnecessary negative environmental impact in connection with use, recirculation and possibly destruction
- Complying with the applicable legislation
- Ensuring ongoing improvement of our environmental management system and associated environmental performance

All subsidiaries in the Guldmann group are covered by the above policy, and we expect that our Partners (suppliers and distributors) live up to this policy.

All Guldmann employees are obliged to immediately inform the management if they become aware of any violation of the environmental policy internally in the organisation or at our Partners.

This considers the economic and technological resources at our disposal and our general financial goals for the company and based on our fundamental values.

7.00 Warranty and service conditions

A. Warranty

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the “Warranty Period”). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

Only for USA

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann’s Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

B. Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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2. GB/US

The user lies on his/her back.

Fold either wide end just at the seam, in the direction of the label. Then fold either side in toward the middle of the sling. Fold the remaining side in toward the middle of the sling, overlapping the first side.

Place the folded end, label side down, onto bed or other care surface. Place hand, palm down on top of folded end, to compress bed/surface and slide beneath patient lumbar area.

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Guldmann Sverige AB
Guldmann Srl
Guldmann Inc.