DK	Planløftebåre 255 kg inkl. løftestykke/løftestropper	2
GB/US	Horizontal lifter 255 kg incl. lifting sheet/lifting straps	9
DE.....	Horizontaler Liegebügel 255 kg einschließlich Liegetuch/Hebegurte.....	17
SE.....	Horisontell planlyft 255 kg inkl. lyftlakan/lyftband	25
NO	Planløftebåre 255 kg inkl. løftestykke/løftestropper	32
FR.....	Lève-personne horizontal, 255 kg avec harnais de levage/sangles de levage	39
IT	Sollevatore orizzontale da 255 kg con telo di sollevamento/cinghie di sollevamento ...	47
NL.....	Horizontaal juk 255 kg incl. draagband/tilband	55

User manual – vers. 102.00

Guldmann™

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Item nos:

28444 Horizontal lifter, foldable – 2 strap version

28445 Horizontal lifter, foldable – with weight adjustment

1.00	Use – Horizontal lifter	9
1.01	Manufacturer	9
1.02	Intended purpose	9
1.03	Area of use	9
1.04	Conditions of use	9
1.05	Important/warnings	9
1.06	Unpacking	10
1.07	Before you can use the product	10
1.08	Using the horizontal lifter	10
2.00	Description of features	11
2.01	Labels	11
2.02	Description of feature	11
2.03	Installing the horizontal lifter in the hoist	11
2.04	Accessories	12
3.00	Transport	12
4.00	Maintenance and storage	12
4.01	Cleaning and maintenance	12
4.02	Storage and use	12
4.03	Owner's daily maintenance obligations	13
4.04	Lifetime	13
4.05	Disposal of the horizontal lifter	13
5.00	Technical specifications	13
6.00	Use – Lifting sheet/lifting straps	13
6.01	Intended use	13
6.02	Area of use	13
6.03	Conditions for use	13
6.04	Important/Precautions	13
6.05	UHF RFID tag	13
6.06	Labels and Marking	13
6.07	Use	14
7.00	Maintenance	14
7.01	Cleaning	14
7.02	The owner's daily maintenance duty	14
7.03	Disposal of slings	14
8.00	Service and lifetime	14
8.01	Safety/service inspections	14
8.02	Lifetime	15
9.00	Technical specifications	15
10.00	EU-Declaration of conformity	15
11.00	Environmental policy statement – V. Guldmann A/S 15	
12.00	Warranty and service conditions	16
A.	Warranty	16
B.	Service or Repair	16
	Placing the Lifting sheet	62
	Placing the Lifting straps	74

1.00 Use – Horizontal lifter**1.01 Manufacturer**

V. Guldmann A/S
 Graham Bells Vej 21-23A
 DK-8200 Aarhus N
 Tel. + 45 8741 3100
 www.guldmann.com

1.02 Intended purpose

The Lifting hanger is intended for lifting and transferring a person with disabilities in a horizontal position.

1.03 Area of use

The Horizontal lifter is suited for use in hospital, nursing homes and in institutions.

1.04 Conditions of use

The Horizontal lifter is designed for use with a ceiling hoist systems and may be used only for horizontal relocation* of people lying on a level, horizontal surface. A level surface can be, for example, a bed, a stretcher or the floor.

** A horizontal relocation is defined as a relocation of a person lying on a horizontal surface where the horizontal position is kept during a vertical, horizontal or combined movement. The surface can be completely level or at an angle suited to the respective situation.*

The use of the horizontal lifter is subject to the following prerequisites:


- The device/working tool is used by trained personnel.
- The assistant pays attention to the well-being of the user/patient when using the product.
- The product may be used only with a hoist and rail system which have been tested and approved according to Guldmann's instructions.

1.05 Important/warnings

Always exercise caution when using the horizontal lifter. This means that assistants should always make sure that they, the user and any other persons are safe so as to avoid injuries or accidents.

- Please read the user manual before you start using the horizontal lifter.
- Do not exceed the horizontal lifter's recommended load of 255 kg (560 lbs).
- Products marked with the lowest carrying capacity determine the maximum load.
- Use the horizontal lifter only to move people "in a horizontal position".
- If a defect appears during use of the horizontal lifter, stop using it and contact the Guldmann Service Team for repairs.
- The user in the horizontal lifter must always be under supervision and may not turn or move in the lengthwise direction when lying in the lifter.
- To avoid unintentional movements of the horizontal lifter, you must position the hoist centrally over the horizontal lifter before you start lifting.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.
- Do not modify this equipment without the authorization of the manufacturer

- When the horizontal lifter is lowered with a user, the lifter may never rest on the user. This could happen in two situations:
- When the user is unloaded onto a bed which is narrower than the lifter 770 mm (30 3/4")
The horizontal lifter could continue lowering by the side of the bed and thereby rest on the user (see figure 1).
- When the user is unloaded onto a bed which is wider than the lifter 770 mm (30 3/4")
The horizontal lifter could rest on the bed which the user lies on and thereby open up so that the central part rests on the user. The horizontal lifter can, however, remain closed during unloading if the lifting sling or the sling straps are firmly tightened. If the straps are loosened in this situation, the lifter will open suddenly and rest on the user. (see figure 2).

1.06 Unpacking

Visual check of the horizontal lifter

If the horizontal lifter is thought to be damaged upon receipt, it must not be used before it is checked and approved by a qualified person or the Guldmann Service Team.

The contents of the box


- Horizontal lifter
- Oval hook
- Lifting sheet, incl. bracing struts
- Wall bracket (storage)
- Manual

1.07 Before you can use the product

The horizontal lifter is a foldable working tool which can be suspended in a hoist and unfolded before use. It can generally be stored in two ways:

- On a wall bracket, hanging down the wall.
- In a specially furnished transport trolley where it can be transported to alternative places of use or to storage/cleaning.

There is a risk of getting your fingers caught during the unfolding of the horizontal lifter.


Set up from a wall bracket (storage)

Attach the snap hook of the hoist to the lifting eye of the horizontal lifter manually and unfold the lifter by raising the hoist.

Be careful not to get your fingers caught during unfolding.

When the horizontal lifter unfolds fully, it is ready for use.


Set up from a transport trolley (accessory)

Position the hoist centrally over the horizontal lifter in the transport trolley. Attach the snap hook of the hoist to the lifting eye of the horizontal lifter manually; the lifter will be raised by the hoist and will unfold by the lifting movement.

Be careful not to get your fingers caught during unfolding.

When the horizontal lifter unfolds fully, it is safe to use.

1.08 Using the horizontal lifter

Attention should be paid to the functional level of the user and the work technique of the assistant.

The user should be lifted only to the extent that she/he is clear of the surface and then moved at this height.

The horizontal lifter is used with a Lifting sheet/lifting straps with 12 straps that are fitted in the strap locks. The individual fastening points must be used in pairs and symmetrically.

Lifting straps are an additional accessory (28448).

For proper fitting of the Lifting sheet/lifting straps, refer to 62 and 74.

Slings made by other manufacturers

Guldmann disclaims any responsibility for defects or accidents which occur when slings made by other producers are used.


In case of doubt concerning the selection or use of a lifting sling, please contact your supplier.

Guldmann disclaims any responsibility for defects or accidents which occur as a result of **incorrect** use of the sling or because of a lack of attention from the carer or user.

2.00 Description of features

2.01 Labels

ID number label

Injury warning (yellow)

2.02 Description of feature

Fig. A

Foldable horizontal lifter with centre of mass adjustment

1. Lifting eye for Snap hook (attachment to hoist)
2. Handle for centre of mass adjustment (horizontal position)
3. Folding (saves room when not in use)
4. Telescopic tubes (save room, length adjustment sheet/user)
5. Strap lock (length adjustment sling/strap)
6. Lifting sheet/strap set
7. Telescopic arms

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


2.03 Installing the horizontal lifter in the hoist

The lifting eye of the horizontal lifter must always be attached securely to the ceiling hoist during preparation/use:

1. GH3 – the lifting strap can be fitted using the oval hook supplied, (see figure 1).
2. GH3 – the lifting strap can also be fitted using a swivel (accessory), (see figure 2).

The lifting hook is fitted by moving the hook opening across the flat part of the hanger on the swivel. When the hook cannot be moved further onto the hanger, twist it into place in an upright position. The horizontal lifter can be used when the lifting strap has been moved up and the hook is in an upright position.

3. GH2 – the snap hook on the lifting strap must always be fitted in such a way that the lock plate is engaged, (see figure 3).

When the straps of the hoist are not long enough for the desired transfer or lift, a Guldmann extension strap can be used (see section 2.03).

Lifting sheet/lifting straps, strap lock

The horizontal lifter is used with a Lifting sheet/lifting straps with 12 straps that must all be securely attached to the strap locks.

Placing the Lifting sheet/lifting straps, see 62 and 74.

Centre of mass adjustment

The horizontal lifter is equipped with a centre of mass adjustment, refer to item on Page 8, it. 2 on the main drawing. You should strive to place the user in the horizontal lifter in such a way that the lifter is level when the centre of mass adjustment is in a neutral position so as to achieve an even load on both lifter and sling.

You can adjust the centre of mass of the horizontal sling by turning one of the hand wheels placed in each end of the center part. By turning the wheel clockwise, you will raise the horizontal lifter in the end where you turn. By turning the wheel counter-clockwise, you will lower the horizontal lifter in the end where you turn.

Lifting up and down

When all straps on the Lifting sheet/lifting straps have been fastened to the strap locks, you can start the lift. Make sure that the centre of mass adjustment is set correctly before continuing with the lift. If the deviation in the angle is significant (+/-15 degrees) or greater than intended, the centre of mass must be adjusted.

If deviations are too large, the horizontal lifter must be lowered and the user must shift position on the sheet.

Movement

For all movements up/down/horizontally, you must make sure that the horizontal lifter/ sling/user does not come into contact with their surroundings in a way that could give rise to dangerous situations, for example, that the horizontal lifter is brought out of balance or impacts, gets caught in or overturns other objects.

2.04 Accessories

Extension strap

When the straps of the hoist are not long enough for the desired movement or lift, you can use a Guldmann extension strap. This can happen, for example, when the hoist system is placed high and you need to be able to lift the user from the floor.

Approved Guldmann extension straps:

	Length
126845-100	10 cm
126845-150	15 cm
126845-300	30 cm
126845-400	40 cm
126845-500	50 cm
126845-600	60 cm
126845-700	70 cm
126845-800	80 cm
126845-900	90 cm
126845-1000	100 cm

Sling set

Lifting straps, part number 28448, comprising
1 head support
6 lifting straps
1 sling pin

Placing the lifting straps, see 74.

The sling set is primarily intended for use at intensive care units and for persons who cannot tolerate being turned over. It is also the case with this sling set that the weight must be distributed among all 12 straps in order to lift the max. load. The individual fastening points must be used in pairs and symmetrically.

Transport trolley

The transport trolley, part no. 28443, is adapted for the transport of the horizontal lifter to alternative locations for use or storage/cleaning.

Swivel

Swivel GH2/GH3, spare part – part no. 284520

When using the swivel in place of the oval hook, the horizontal lifter can be turned through 360 degrees during use without the lifting strap becoming twisted.

3.00 Transport

Guldmann recommends that the horizontal lifter is transported or stored in its original packaging. Alternatively, the horizontal lifter can be transported or stored in a transport trolley, which is an accessory to the horizontal lifter.

4.00 Maintenance and storage

4.01 Cleaning and maintenance

Clean the horizontal lifter with a damp cloth and ordinary detergent. Do not use strong chemicals or an autoclave.

The horizontal lifter does not need any additional maintenance apart from ordinary cleaning.

4.02 Storage and use

The horizontal lifter is designed for indoor use.
Temperature +5°C to +40°C
Humidity <95% (not condensing)

For prolonged storage, the horizontal lifter should be stored in a dry room, preferably suspended from the wall or in the transport trolley.


Short-term storage in between uses of the horizontal lifter is described below.

Wall bracket (storage)

For storage of the horizontal lifter after use, you can use a wall bracket mounted to a solid vertical wall. The horizontal lifter is then suspended down the wall. The wall bracket must be placed as close to the floor as possible in order to ensure a good working height – the horizontal lifter must be able to unfold completely and hang free.

Storage on wall bracket

Retract the telescopic tubes fully to minimise size. Then place the two central strap locks in the middle of the horizontal lifter. Place one of the telescopic arms of the horizontal lifter into the wall bracket by guiding the "hanging" horizontal lifter into the wall bracket manually. Then lower the horizontal lifter until it hangs freely down the wall and the snap hook can be released from the lifter.


Transport trolley

If you want to use the hoist at several locations or store it centrally, the horizontal lifter can be placed and transported folded in a transport trolley.

Storage in transport trolley

Fully retract the telescopic arms to minimise size and place the two central strap locks in the middle of the horizontal lifter. Place the trolley centrally under the horizontal lifter with the arms spread out so the horizontal lifter can be placed into the slot manually while being lowered. The horizontal lifter is then folded in the trolley by means of the lowering movement.

Do not touch the central part and the hinges of the horizontal lifter during folding, as your fingers can be caught in the moving parts.

The horizontal lifter is safe to handle when it is fully folded.

4.03 Owner's daily maintenance obligations

Check the sling and the sling straps for wear and tear and damages before using them. Do not use the sheet if it is damaged or defective.

Do not use the horizontal lifter if the strap locks or other mechanical parts are damaged, loose, defective or missing. Instead, contact your supplier for a service visit.

Make sure that the screws on the horizontal lifter are tightened and that there is no abnormal backlash. Call for a service technician if there is any backlash.

Pay special attention to:

- Screws at the end cover of the central part.
- Screws in the strap locks + the strap lock bracket.

Replacement of components

Components may be replaced by a qualified service technician or the Guldmann Service Team.

4.04 Lifetime

The horizontal lifter has an expected useful life of 10 years, when subjected to ongoing maintenance.

4.05 Disposal of the horizontal lifter

Comply with the local and national regulations on environmentally correct recycling. All parts can be separated for recycling.

5.00 Technical specifications

Lifting capacity max. 255 kg
Weight 26 kg

Dimensions

Horizontal lifter, unfolded (W x H x L) 1800 x 970 x 490 mm
Horizontal lifter, parked (W x H x L) 1100 x 230 x 1215 mm

Materials

Horizontal lifter. Aluminium, anodised

6.00 Use – Lifting sheet/lifting straps

6.01 Intended use

The Lifting sheet/lifting straps is intended for lifting or supporting a person or body parts of a person

6.02 Area of use

The Lifting sheet/lifting straps is suited for use in hospitals and in nursing homes.

6.03 Conditions for use

The Lifting sheet/lifting straps is designed for use with a horizontal lifting stretcher in combination with a ceiling hoist system. It may be used only for a horizontal relocation of a person lying on a level, horizontal surface as for instance a bed, stretcher or the floor.

The Lifting sheet/lifting straps is used for horizontal relocation procedures involving people with reduced mobility.

The use of the Lifting sheet/lifting straps is subject to the following:

- The Lifting sheet/lifting straps is used by trained staff or persons who have been instructed in the use of the Lifting sheet/lifting straps in question.
- The maximum lifting capacity of 255 kg must never be exceeded.
- The Lifting sheet/lifting straps is to be used for horizontal transfer of a user lying on a level, horizontal surface.
- The Lifting sheet/lifting straps must be used with a Guldmann horizontal lifter.
- The Lifting sheet/lifting straps must be secured with all 12 straps fitted to the strap locks on the horizontal lifter.
- The carer must pay attention to the user's well-being when using the Lifting sheet/lifting straps.

6.04 Important/Precautions

- Read the instructions carefully before using the Lifting sheet/lifting straps.
- The Lifting sheet/lifting straps's maximum load must never be exceeded.
- The Lifting sheet/lifting straps may only be used to transfer a user who is in a horizontal or angled position when using a horizontal lifter.
- Before a Lifting sheet/lifting straps is used, it must be examined according to point 2.02.
- Possible repairs must only be made by the manufacturer.

6.05 UHF RFID tag

This product features a passive UHF RFID tag integrated in the product label. The RFID tag can be used for asset management & tracking purposes.

The RFID tag is readable with equipment compliant to EPC global UHF Class I Gen ISO 18000-63.

6.06 Labels and Marking

CE marking

Medical Device Class I in accordance with EU MDR Regulation



Read the manual before use


UK Responsible Person

European Device Solutions Ltd. 15 Coanwood Drive, Whitley Bay, Tyne & Wear, NE25 9GB, United Kingdom.
Email: info@europeandevicesolutions.co.uk
Tel: +44-754-559-5464


Example of serial number label Inspection label

Product label

Size label

Tumble-drying at low temperature

Do not iron

7.02 The owner's daily maintenance duty

Check the lifting sling for wear and damage before use according to the following checklist which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Sling inspection checklist

Before using a Guldmann sling / accessory check the following:

Is the sling clean?

Follow facility specific infection control procedure.

Is the sling's label present, legible and complete?

Missing, illegible or incomplete sling label(s) could make identification of appropriate size of the sling, function of sling, and or weight limit capacity of the sling impossible.

Are the lifting straps and stitches intact?

- Look for broken or worn stitches
- Look for knots in straps
- Look for tears or fraying of straps
- Look for snags or punctures or holes
- Look for any particles in fabric or straps

Is the fabric intact?

- Look for abnormal wear patterns, excessive wear, abrasive evidence
- Look for cuts or frayed fabric
- Look for unusual or significant discoloration
- Look for snags, punctures, tears, holes
- Look for frayed or insecure seams
- Look for any acid / caustic / thermal burns
- Look for changes in material consistency, e.g. increased stiffness
- Look for any imbedded particles

Has the shape of the sling been altered, made shorter or longer in relation to the original size using knots, needles, tape or other methods?

Conclusion

If the sling suffers from one or more of the above mentioned conditions then it must be taken out of service regardless of the weight of the person to be lifted.

6.07 Use

If there is any doubt about the selection or use of a Lifting sheet/lifting straps, please contact your supplier.

Important!

Plan the transfer and do not leave the user unsupervised in the Lifting sheet/lifting straps. Before lifting, ensure that the patient cannot become stuck and that the Lifting sheet/lifting straps is not caught on the bed, wheelchair or on other surrounding items. The patient's head, arms, hands and feet must not be in danger of becoming trapped. Be careful with any tubes or wires that may be attached to the patient and/or nearby equipment. Check that the hand control and hand control cable are free of the horizontal lifter, user and other objects before the lifter is moved up or down.

Be careful when attaching the Lifting sheet/lifting straps. Ensure that the straps have been correctly fitted in the strap locks on the horizontal lifter before the user is lifted, refer to description under fitting of Lifting sheet/lifting straps. Before lifting, ensure that the straps remain in the correct position.

Guldmann shall not be liable for faults or accidents due to incorrect use of the Lifting sheet/lifting straps, or for reasons of inadequate attention on the part of the carer or user. If the Lifting sheet/lifting straps is used in combination with products that are not manufactured by Guldmann, a risk assessment must be made by qualified staff.

Placing the Lifting sheet/lifting straps, see 62 and 74.

7.00 Maintenance

7.01 Cleaning

Normal washing at the indicated temperature

Do not use bleaching agent

7.03 Disposal of slings

Slings are disposed of by incineration. By proper incineration polyester will be degraded to carbon dioxide and water.

8.00 Service and lifetime

8.01 Safety/service inspections

In accordance with international standard EN/ISO 10535 "Hoist for the transfer of disabled persons – Requirements and test methods" an inspection **must** be performed every 6-month according to the following instructions, which is not intended to represent all potential inspection steps. Potential damage may vary. Judgment of inspector/site prevails.

Safe Operating Practices with Slings

Considerations for damaged or defective slings and taking them out of service:

Withdraw the sling from service if one or more of the following conditions are present:

- chemical or caustic burns
- melting or charring of any part of the sling
- snags, punctures, tears or cuts
- broken or worn stitches
- missing, illegible or incomplete sling tag
- knots in any part of the sling
- abrasion
- other visible damage that causes doubt as to the strength of the sling

Sling inspection is done for the protection of the user, the caregiver, and the overall hospital site safety. A sling inspection system has additional benefit. Systematic sling inspection will assist in the identification of damage trends, potentially leading to cost effective suggestions and results. The inspection process can also help to identify inventory duplicity in certain sling types and sizes.

Sling inspection system

Development of a specific procedure and program for the inspection of slings at your facility is your best safeguard. Consider employing a three part system of inspection. Slings that are removed from service and are not capable of repair should be disposed of so they are unfit for any future use and can not find a way back into active inventory.

1) Initial

This level of inspection is done at the time that the sling is received into your facility. The inspector should ensure that no damage has occurred during transit, and also verify that the sling work load limits match those contained in the manufacturer's catalogue. If your facility documents the sling inspection process through written inspection records, the paper trail should begin at this stage.

2) Frequent

The frequent level of inspection should be done by the sling user before each use. The sling should be examined and removed from service if damage is detected. The sling user should also determine that the sling is proper for the user conditions, care task required and the required weight capacity.

3) Periodic

Your facility might want to consider implementing a program for a periodic level of inspection at regular intervals. The interval should be based upon the frequency of use, severity of the service cycle and information derived through the inspection process. Recommendations to prevent damage and enhance service life could be made by staff that perform the periodic inspections. If written inspection records are maintained, they should always reference the unique sling identification number, and be updated to record the condition of the sling. Not intended to represent all potential inspection steps or all potential aspects of product management program. Judgment of inspector/site prevails.

Sling inspection technique

The sling inspection procedure should be thorough, systematic and consistent; both visual and "hands on" inspection techniques are recommended. Certain forms of damage are far more discernable through hands-on inspection, than by visual inspection. For example, fabric stiffness, crushed webbing, as well as, thinning fabric can be identified through tactile inspection. Visual inspection alone may not reveal all forms of sling damage. Once signs of damage have been identified, do not downgrade the work load limit of the sling, with the intent of continuing to use it, but at limited capacity or frequency. This is sometimes done to get more service life out of a damaged sling. The operating rule and standard should be: intact = use; damage = do not use.

Consider the practice of documenting sling inspections through written inspection records.

The documentation should include information such as: the name of manufacturer, the sling stock number, width and length, the

unique sling identification number (important in differentiating similar slings), as well as the condition of the sling. Other important information might also include the date it was received or put into use at your facility and any special features (if applicable). A beneficial outcome of an inspection program would be the realization of repetitive forms of damage and the analysis that would lead to specific recommendations.

Sample visual examples of synthetic sling damage ^{x)}

Chemical/caustic burns

Broken stitching

Crushed / Frayed webbing

Knots

Melting / Charring

^{x)} sample visual images not intended to represent all types of potential damage

8.02 Lifetime

The life of the sling is individual and depends on how it is used, washed etc. Before use the sling must be examined according to description in section 2.02 and if it does not meet the inspection requirements, it must be discarded if necessary.

9.00 Technical specifications

Lifting capacity, SWL 255 kg (560 lbs)
Material Polyester

10.00 EU-Declaration of conformity

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

11.00 Environmental policy statement – V. Guldmann A/S

Guldmann is continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

It is Guldmann's goal to:

- Comply with the current environmental legislation (e.g. WEEE and REACH directives)
- Ensure that we, at the widest possible range, use RoHS compliant materials and components

- Ensure that our products do not have an unnecessary negative impact on the environment regarding use, recirculation or disposal
- Ensure that our products contribute to a positive working environment in the places they are utilised

Inspections are made annually by the Department for Nature and Environment from the Municipality of Aarhus using the Danish Environmental Protection Act, section 42 as a reference.

12.00 Warranty and service conditions

A. Warranty

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

The warranty does not cover any part of the equipment that has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment that has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann certified technician. Any parts or components repaired or replaced by a Guldmann certified technician will be guaranteed for the remainder of the Warranty Period.

Only for USA

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann Certified Technician. A Guldmann Certified Technician is a technician who has successfully completed Guldmann Service Training, and who holds a valid Service Training Certificate from Guldmann, and is in possession of a valid password to access Guldmann's Service and Information Console (SIC). A Guldmann Service Training Certificate and SIC password are valid for three years (only USA) from the date the technician is first certified. Thereafter, the technician must undergo re-certification training to obtain a new valid certificate and password. Any parts or components repaired or replaced by a Guldmann Certified Technician will be guaranteed for the remainder of the Warranty Period. In the event the warranty is rendered null and void, the purchaser shall indemnify and hold Guldmann harmless of and from any and all claims or liability arising as a result of equipment malfunction or misuse.

B. Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

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1. GB/US

The lifting eye of the horizontal lifter must always be attached securely to the ceiling hoist during preparation/use, either using the oval hook, swivel (accessory) or using the GH2 hoist snap hook, (see section 2.02).

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
2. GB/US

Lifting sheet, strap lock

The lifting sheet has 12 straps that must all be securely attached to the strap locks.

Make sure that you turn the lifting sheet in such a way that the labelling of the lifting sheet faces away from the user/patient.

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3 GB/US

The straps are infinitely adjustable, however the distance from the seam on the end of the strap to the lock must always be at least 30 mm.

The straps must be fastened to the strap locks by opening the locks, passing the straps through and letting go of the lock. The strap lock will close automatically. Ensure that the strap is locked correctly. The straps must always be fastened from the inside and towards the carer.

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4. GB/US

Fold the lifting sheet into a Z shape against the centre marking of the lifting sheet so that it turns away from the patient.

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5. GB/US

Place the lifting sheet alongside the patient on the surface.

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Vend patienten i sideleje og placer sejlet langs patientens ryg.

6. GB/US

Turn the patient into a lateral position and place the lifting sheet along the patient's back.

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7. GB/US

Turn the patient back onto the lifting sheet. The lifting sheet can now be unfolded. If necessary, you can turn the patient onto their opposite side before unfolding the lifting sheet.

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Kontroller at patienten ligger midt på sejlet.

8. GB/US

Check that the patient is lying in the middle of the lifting sheet.

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9. GB/US

If necessary, insert the side stiffeners in the two pockets.

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10. GB/US

The side stiffeners must be connected in the middle by screwing them into the accompanying threaded piece.

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11. GB/US

Fit the lifting strap onto the horizontal lifter.

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16. GB/US

When the patient is placed on the bed again, the lifting sheet can be removed by taking the straps out of the strap locks on the horizontal lifter.

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17. GB/US

The stiffeners can be removed from the pocket at the main end or they can be split in the middle by unscrewing one of the stiffeners from the threaded piece and then pulling the two parts out of the pocket in the middle of the sling.

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
18. GB/US

The stiffeners can now be pulled out of the pocket.

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
19. GB/US

For cleaning

Unscrew the threaded pieces from the stiffeners.

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GB/US

Placing the Lifting straps
(additional accessory, 28448)

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1. GB/US

Insert the provided plastic stick into the pocket on the lifting strap.

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2. GB/US

Rotate the lifting strap so that the plastic stick faces downward towards the mattress. Now, slide the lifting strap and plastic stick underneath the user. After placing the lifting strap, remove the plastic stick. Repeat the same process with the remaining straps.

Ensure that the straps are facing in a way that the markings on the straps are away from the user/patient.

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Guldmann GmbH
Guldmann Sverige AB
Guldmann Srl
Guldmann Inc.