



Guldmann™

## Guldmann Limb Sling – Quiz

**1. Which care tasks would you select the Guldmann Limb Sling to help you perform:**

- A. Assessment of patients skin on limbs, feet, heels upon admission and as indicated
- B. Elevate legs to place pillows to take pressure from patient's heels
- C. Wound dressing changes or wrapping of limbs
- D. Position and support of both limbs to assist with perineal care tasks
- E. All of the above

**Answer E**

**2. What is the pocket of the Limb Sling used for?**

- A. To hold Instructions for sling use
- B. To place hand, so sling can be guided beneath leg with the recommended placement technique
- C. To tuck folded end of the sling into
- D. B and C

**Answer B**

**3. How do you know when a Guldmann Limb Sling should be taken out of service?**

- A. One of the loops is missing from the Limb Sling
- B. After it is has been in use for 6 months
- C. The label is illegible on the Limb sling
- D. A and C

**Answer D**

**4. You are using a Limb sling and becomes soiled during a care task, what do you do with the Guldmann Limb sling when the care task is complete?**

- A. Single patient use slings can be discarded
- B. Reusable/washable Limb slings are put into the appropriate linen hamper
- C. Always discard all slings after use
- D. Use sling for the next patient
- E. A and B

**Answer: E**

**5. When do you explain use for and purpose of the Limb sling to your patient?**

- A. Only when you have time to give an explanation
- B. Only when there are family members in the room
- C. Only for patients that may complain about the use of slings
- D. None of the above

**Answer: D**

- 6. What should you do if you need to use the Limb sling to re-do dressings on a patient's leg but you do not remember how to use it properly?**
- A. Ask your strongest co-worker to help manually lift and then hold up the patient's leg
  - B. Ask your unit peer leader or co-worker to review how to apply and use the Limb sling with you for the task
  - C. Leave the dressing for someone that knows how to use the Limb sling
  - D. Do a quick review of the instructions and video for this sling which are available on the available learning platform or Guldman Academy
  - E. B and D

**Answer: E**

- 7. I can use the Limb sling to lift and keep the patient's leg(s) elevated for an extended period of time**
- A. True
  - B. False

**Answer: B False**

- 8. How many Limb slings should I use to position limbs for perineal access?**
- A. One Limb Sling per leg that needs support and holding
  - B. Two Limb Slings per leg that needs support and holding
  - C. You can place both legs in one Limb sling for perineal access
  - D. It is inappropriate to use Limb Slings for perineal access

**Answer: A**

- 9. Where is the appropriate anatomical location to place the limb sling for the patient's leg?**
- A. Under the thigh to lift the whole leg
  - B. Under the calf to lift the leg while allowing knee and hip flexion
  - C. Above the posterior ankle to lift the leg with no knee flexion and some hip flexion
  - D. B and C

**Answer: D**

- 10. What should I consider if my patient's leg has a much larger proportion than the Limb sling?**
- A. Consider using 2 Limbs Slings adjacent to each other to make the support surface wider for patient's leg
  - B. Consider using a Guldman Multi Support sling as an alternative to the Limb Sling
  - C. Pad the Limb sling with a soft material to increase the comfort and surface area for the patient's leg
  - D. All of the above

**Answer: D**