



# Guldmann™

## Instructions Basic High



### Applying in Bed



Place sling on user's shoulders and or the bed behind the user's head & shoulders. To ensure the sling is centered, the center stitching marks and the pocket should follow the user's spine.




Have user lean forward or assist as needed. Put your hand in the sling pocket and slide sling to follow the back all the way down to the bend of the bed.



Hold your hand between the user's hip and in the inside of the sling at the back support and press down. At the same time pull the leg support down toward the knee with the other hand.



Place the leg support under the thigh. Repeat on the other side/leg. Cross the leg straps (through middle) before attaching to the hanger.



To lift user evenly, match the colors on the loops of each strap, such as grey to grey or green to green. To modify user's position or comfort, differing loop colors can be used - consult with your Guldmann representative for assistance.



Do not lift higher than necessary. To put user into a chair, just before and as one lowers the user, apply pressure to the user/sling against the chair-back in order to get the user as far back into the chair seat as possible. Make sure to use the correct lifting hanger.

On/off in wheelchair

