



Guldmann™

Instructions Net Repositioning



Hanger parallel orientation



Hanger perpendicular orientation

Using Net Repositioning Sling to Lift Full Body



Place the sling, label side down, and at the top, near patient's head. The patient's head should be supported by the sling. The lower portion of the sling should support the patient at a point between their ankle and feet depending upon their height.

If using the standard hanger bar, there are two options for position of hanger bar. When the hanger bar is perpendicular (or across the body), the pt might have more comfort across their chest or trunk. When the hanger bar is parallel to the pt (lengthwise with body), taller pt's are better supported throughout their length.



Step #1:
With the hanger parallel or perpendicular, attach straps # 2 & 3 on each side, using the same color loops. Raise the hanger to reduce slack in lifting straps.



Step #2:
Now take strap # 1, then # 4, and stretch out to determine which color loop(s) is best to select for the user's desired position. Check that all straps/loops are on securely.



Step #3:
Raise the user to the desired clearance off the bed surface. Move user toward head of bed to reposition upward; or move user off of bed surface to other desired surface such as stretcher or litter.

Using the Repositioning Sling to Logroll to Sidelying

With the hanger parallel to user, attach straps # 2 & 3 on the side OPPOSITE of where user will be turned; use the same color loops. (If needed for support, attach straps # 1 & 4).

Take care to support the user when rolling, adjusting position of legs, arms, and head if user can not assist. Use siderails, and stand in the safest position to directly observe and guide. Raise hanger to achieve desired sidelying position. Place pillows behind sling/user as desired; then lower the hanger and remove straps. (Disregard pillow placement if purpose for rolling is for access).



Using the Repositioning Sling to Lift Head/Trunk

To create separation between mattress and head/trunk, such as x-ray cartridge placement for portable x-ray - Attach strap #s 1 & 2 on each side, then raise hanger to desired height. Caution - raising hanger too high will displace sling.



Notes:

- Keep pt's comfort in mind
- Place sling under the patient, so user's head trunk, pelvis and legs are supported. Align the landmark stitching at top and bottom with user's spine
- Position the patient for turning by crossing their legs and bringing one arm over the other
- Raise the siderails on the bed for safety and place a pillow between the patient and the siderail; the patient will be very close to the rail