
DK	Elastisk Line	2
GB/US	Elastic Line	3
DE	Elastikseil	5
SE	Elastisk lina	6
NO	Elastisk line	8
FR	Ligne élastique	9
IT	Banda elastica	11
NL	Elastische lijn	12
ES	Línea elástica	14

User manual – vers. 101.0

DK... ELASTISK LINE

Vers. 101.0

Varenr.:
561629 Elastik line

1.00 Formål og anvendelse

1.01 Producent

V. Guldmann A/S
Graham Bells Vej 21-23A
DK-8200 Aarhus N
Tel. + 45 8741 3100
www.guldmann.com

1.02 Erklærede formål

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1.02 Intended purpose

The elastic line is used for mobilisation and training of users. The purpose is to:

- Maintain and train the user's functional level
- Minimise the risks typically associated with long-term immobility.

1.03 Area of use

The elastic line can be used at hospitals, assisted living facilities, physical therapy centres and in private homes for the physical therapy of people with reduced functionality.

1.04 Conditions of use

When the line is used, it is important that:

- The line is used only by qualified personnel or individuals who have been instructed in how to use the elastic line.
- The person who is training is able to participate in the process.
- The carer pays attention to the user's well-being during the training.
- The line is not used for lifting and handling people.
- The line is used together with Guldmann lifting modules

1.05 Precautions

- Do not exceed the maximum load described in section 1.06.
- Before use, the elastic line should be examined for defects and wear such as tearing, seizing, wear or ripping. If it does not meet the inspection requirements, it must be discarded immediately.
- Do not modify this equipment without authorisation from the manufacturer
- Any serious incident that occurs in connection with the use of this product must be reported to the manufacturer and the local authority.

Example of a defective line

1.06 Use of elastic line

The elastic line can be used in several ways:

Used directly on the hanger.

Used with connecting straps (Art. 29356).

2.00 Environmental conditions

Operation

The products operational environment:

- Operating temperatures between +10°C and +35°C
- Relative air humidity between 30 and 70%
- Air pressure between 700 and 1060 hPa

Beside temperature, the same environmental conditions apply for transportation and storage.

- Transport and storage temperatures between -10°C and +40°C

3.00 Cleaning

Clean the line with a damp cloth and ordinary detergent. Do not clean using chemicals or autoclaving.

4.00 Lifetime

The expected useful life of the line is 5 years, depending, however, on use pattern, cleaning, etc. Before use, the line should be inspected as described in section 1.05 and be discarded if it does not meet the inspection requirements.

5.00 Technical specifications

Maximum load on the elastic line 255 kg
 Material (outer) Polyester
 Material (inner) Latex rubber

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Guldmann Sverige AB
Guldmann Srl
Guldmann Inc.