



**V. Guldmann A/S**  
 Tel. +45 8741 3100  
 Fax +45 8741 3131  
 info@guldmann.com  
 www.guldmann.com

## Bath chair & bath stretcher Technical specifications

### Functions

Lifting capacity, max. : 255 kg

### Weight and materials

Complete chair : 25,5 kg

Frame : Rustproof steel  
 Bearings : Rustproof steel  
 Bolts : Rustproof steel  
 Seat : Moulded polyester  
 Footrest : Painted, moulded aluminium  
 Footplate : Moulded plastic  
 Wheels : Rubber with brake, rustproof steel  
 Safety belt : Nylon

### Measurements

A. : 900 mm  
 B. : 500 mm  
 C. : 500 mm  
 D. : 125 mm  
 E. : 535 mm  
 F. : 1600 mm

Item Number/Model : 14600: Bath chair

### Classification

according to ISO 9999 : 12 36 12

### Functions

Lifting capacity, max. : 255 kg

### Weight and materials

Complete stretcher : 54,4 kg  
 Bench : 25,5 kg  
 Base : 25,5 kg

Frame : Rustproof steel  
 Bearings : Rustproof steel  
 Bolts : Rustproof steel  
 Cover for stretcher : PVC cover in nylon  
 Wheels : Rubber with brake, rustproof steel  
 Safety belt : Nylon

### Measurements

A. : 1800 mm  
 B. : 650 mm  
 C. : 125 mm  
 D. : 890 mm  
 E. : 1970 mm  
 F. : 360 mm  
 G. : 1415 mm


Angularity : 3°, 13°, 19° and 23°

Item Number/Model : 14610: Bath stretcher


### Classification

according to ISO 9999 : 12 36 12

### BATH CHAIR



### BATH STRETCHER



Time to care