

**Material**

Panel	: Graphite fibre/Glass fibre
Joint	: Nylon
Handle	: Nylon
Corner protector	: Rubber
Walk surface	: Rubber tip paint
Paint	: Urethane

**Colour**

Walk surface	: Grey
Center Line	: Blue
Side Portion	: Yellow
Rear	: Grey
Velcro tape	: Yellow

**Durability**

Max load	: 300 kg
----------	----------

**Accessory**

Item no.: 30100-000	: Stepless wheelchair carrying bag for 70 cm ramp
---------------------	---

**Marking**

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

**Stepless Lite ramps**

**Item no. 30100-070**

Length	: 700 mm
Width, inner	: 730 mm
outer	: 730 mm
folded	: 365 mm
Height of edge	: none
Thickness (Fold)	: 38 mm
Max height difference	: 140 mm (1:5)
Weight	: 3,5 kg

**Lite ramps**

**Item no. 30100-085**

Length	: 850 mm
Width, inner	: 780 mm (390 x 2)
outer	: 780 mm
folded	: 390 mm
Height of edge	: none
Thickness (Fold)	: 38 mm
Max height difference	: 170 mm (1:5)
Weight	: 4 kg

**Lite ramps**

**Item no. 30100-125**

Length	: 1250 mm
Width, inner	: 760 mm
outer	: 780 mm
folded	: 390 mm
Height of edge	: 25 mm
Thickness (Fold)	: 106 mm
Max height difference	: 250 mm (1:5)
Weight	: 6 kg

**Lite ramps**

**Item no. 30100-165**

Length	: 1650 mm
Width, inner	: 760 mm
outer	: 780 mm
folded	: 390 mm
Height of edge	: 25 mm
Thickness (Fold)	: 106 mm
Max height difference	: 330 mm (1:5)
Weight	: 7,5 kg

**Lite ramps**

**Item no. 30100-205**

Length	: 2050 mm
Width, inner	: 760 mm
outer	: 780 mm
folded	: 390 mm
Height of edge	: 25 mm
Thickness (Fold)	: 106 mm
Max height difference	: 410 mm (1:5)
Weight	: 9,5 kg