



Active Micro Plus			Active Trainer				Gait Trainer		
Size	Item no.	Item no.	Size	Item no.	Item no.	Size	Item no.	Item no.	
Kids 4-6	281001		Kids 4-6	283001		S	283241		
Kids 6-10	281011		Kids 6-10	283011		M	283251		
Kids 10-14	281021		Kids 10-14	283021		L	283261		
XS	281531	281535	XS	283031		XL	283271		
S	281541	281545	S	283041	283545	2XL	283281		
M	281551	281555	M	283051	283555	<b>Gait Trainer Bariatric</b>			
L	281561	281565	L	283061	283565	One size	283100		
XL	281571	281575	XL	283071	283575				
2XL	281581	281585	2XL	283081	283585				
<b>Extra support belt</b>			<b>Active Vest Kids</b>						
XS / S	281533		Kids 2-4	2831000					
M	281553		Kids 4-6	283101					
L	281563		Kids 6-10	283111					
XL	281573		Kids 10-16	283121					
2XL	281583								
3XL	281593								

Active Micro			Basic High			
Size	Item no.	Item no.	Size	Item no.	Item no.	Item no.
XS	284031		Kids 4-6	272001	272103	
S	284041		Kids 6-10	272011	272113	283615
M	284051		Kids 10-14	272021	272123	
L	284061		XS	272031	272133	283635
XL	284071		S	272041	272143	283645
2XL	284081		M	272051	272153	283655
			L	272061	272163	283665
			XL	272071	272173	283675
			2XL	272081	272183	283685
			<b>Basic High Bariatric</b>			
			L / XL			283672
			XL / 2XL			283682
			2XL / 3XL			283692
			3XL	2720913		
			4XL	2720914		
			5XL	2720916		

Basic Basic				Basic Low			Basic Shell	
Size	Item no.	Item no.	Item no.	Size	Item no.	Item no.	Size	Item no.
Kids 4-6	270001	270103		XS	271031	271133	XS	275032
Kids 6-10	270011	270113		S	271041	271143	S	275042
Kids 10-14	270021	270123		M	271051	271153	M	275052
XS	270031	270133		L	271061	271163	L	275062
S	270041	270143		XL	271071	271173	XL	275072
M	270051	270153		2XL	271081	271183		
L	270061	270163						
XL	270071	270173						
2XL	270081	270183						
3XL	270091							

Basic Comfort High				Custom Amputee		Sit-On Comfort High			
Size	Item no.	Item no.	Item no.	Size	Item no.	Size	Item no.	Item no.	Item no.
XS	277031	277032	277035	XS	290031	Kids 4-6	294001		
S	277041	277042	277045	S	290041	Kids 6-10	294011		
M	277051	277052	277055	M	290051	XS	294031		
L	277061	277062	277065	L	290061	S	294041		
XL	277071	277072	277075	XL	290071	M	294051		
2XL	277081	277082	277085			L	294061		
						XL	294071		
						2XL	294081		
						3XL	294091		
						<b>Sit-On Comfort Hygiene</b>			
						XS	294131		
						S	294141		
						M	294151		
						L	294161		
						XL	294171		

Sit-On Comfort			Sit-On II			Sit-On High II		
Size	Item no.	Item no.	Size	Item no.	Item no.	Size	Item no.	Item no.
Kids 4-6	293001		Kids 4-6	297001		Kids 4-6	298001	
Kids 6-10	293011		Kids 6-10	297011		Kids 6-10	298011	
XS	293031		XS	297031		XS	298031	
S	293041		S	297041		S	298041	
M	293051		M	297051		M	298051	
L	293061		L	297061		L	298061	
XL	293071		XL	297071		XL	298071	
						2XL	298081	
						3XL	298091	

## Guidelines for choosing the size of a sling

Max lifting capacity for all Guldmann seated and standing slings: 255 kg (560 lbs)

Exceptions:  
Basic Basic 3XL, Basic High Bariatric 3XL, 4XL, 5XL and Gait Trainer Bariatric 500 kg (1100 lbs)

Active Vest Kids max lifting capacity: 100 kg (220 lbs)


Disposable slings made of polyester or polypropylene 205 kg (450 lbs)

Disposable slings made of polyester and cotton 350 kg (770 lbs)

W2 = Chest measurements

x) Guideline for choosing Active Trainer and Gait Trainer

xx) Guideline for choosing Active Vest Kids



Size / cm	W1	W2 <sup>x)</sup>	W2 <sup>xx)</sup>	H	Size / inches	W1	W2 <sup>x)</sup>	W2 <sup>xx)</sup>	H
Kids 2-4			46 - 54		Kids 2-4			18 1/8 - 21 1/8	
Kids 4-6	25 - 30	55 - 90	50 - 65	45 - 60	Kids 4-6	9 7/8 - 11 5/8	21 1/8 - 35 3/8	19 1/8 - 25 1/8	17 1/8 - 23 1/8
Kids 6-10	28 - 33	60 - 95	51 - 72	55 - 70	Kids 6-10	11 - 13 5/8	25 1/2 - 41	20 - 28 5/8	21 1/8 - 30 3/8
Kids 10-14	31 - 36	65 - 100		68 - 80	Kids 10-14	13 - 15	27 - 42 1/2		26 3/8 - 31 1/8
Kids 10-16			56 - 75		Kids 10-16			22 - 29 1/8	
XS	33 - 37	70 - 105		75 - 82	XS	13 - 14 5/8	28 1/8 - 44		29 1/8 - 32 1/8
S	36 - 40	75 - 110		80 - 87	S	14 1/8 - 15 5/8	30 - 45 1/8		31 1/8 - 34 1/8
M	39 - 43	80 - 120		85 - 92	M	15 1/8 - 16 5/8	31 1/8 - 47 1/8		33 1/8 - 36 1/8
L	42 - 45	90 - 130		90 - 97	L	16 1/8 - 17 5/8	35 1/8 - 51		35 3/8 - 38 3/8
XL	45 - 49	100 - 140		95 - 102	XL	17 1/8 - 19 1/8	39 1/8 - 55		37 3/8 - 40 1/8
2XL	48 - 54	110 - 150		95 - 102	2XL	18 1/8 - 20 1/8	43 1/8 - 59		37 3/8 - 40 1/8
3XL	54 - 58			95 - 102	3XL	21 1/8 - 22 1/8			37 3/8 - 40 1/8
4XL	57 - 63			95 - 102	4XL	22 1/8 - 24 1/8			37 3/8 - 40 1/8
5XL	62 - 73			95 - 102	5XL	24 1/8 - 28 1/8			37 3/8 - 40 1/8

Time to care

Guldmann™