



V. Guldmann A/S
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com

GLS5.2 Active lifter Technical specifications

Functions

Lifting capacity, max : 155/205 kg

Operation

Lift : Electric
Width adjustment : Electric
Knee support : Manual
Lifting boom : Manual

Weight

Totally : 54 kg
Chassis without foot plate
& knee pad : 23 kg
Mast and lifting boom
incl. control box & battery : 20 kg
Hand control : 0.25 kg

Measurements

A : 25 mm
B : 90 mm
C : 1770 mm
D : 1237 mm
E : 575 mm
F : 1210 mm
G : 460 mm
H min/max : 850/1740 mm
I min/max : 350/450 mm
J min/max : 560/1190 mm
K min/max : 670/1310 mm
L min/max : 230/320 mm

Turning radius

Turning radius : 1360 mm

Safety Features

Emergency stop : Yes
Emergency lowering : Yes, electronic. Mechanical on request
Protection against trapping : Yes

Electrical parts

On/off : Automatically
Integrated Power supply for charging
Input : 120-240Vac, 50-60 Hz, 30 VA
Battery, replaceable : 25.2 V / 2.9 Ah
Charging time : Max. 4 hours
Consumption/power of actuator : 24V, max 8A
Duty Cycle : Max 10%, max. 2 min on, 18 min off.
Battery protection : Switches off automatically at low voltage. Diode light up to indicate critical voltage.

Class of tightness

Mobile hoist : IP 30
Hand control : IPX 4
Integrated Power supply : IPX 6
Battery : IPX 6
Actuator : IPX 4

Accessories/

Miscellaneous : External charger

Labelling



The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

562200 : 155 kg
562201 : 205 kg

Classified

Acc. to ISO 9999 : 12 36 03



Time to care