



Smart Glide

Making daily tasks easy, safe and efficient

Smart Glide is an all-round sliding sheet used for reducing friction between the user and the surface. It can be used for different care tasks providing a more efficient and safe working process.

The sliding sheet is a great supplement to our sling range and eases the process of positioning the sling before lifting the patient either from wheelchair or bed. This ensures that the movement is safer and more comfortable for both patient and care staff.

Smart Glide is made from nylon material with a silicone coating that minimizes friction and noise.

The fabric is OEKO-TEX 100 certified and thus tested for hazardous substances harmful to both health and environment.

For disinfection and cleaning Smart Glide can either be wiped off with disinfectant or laundered.

The Smart Glide fabric has a max weight capacity of 150 kg.



- Supports easy use of slings
- No noise during usage
- Cut only the size you need
- Designed with rip-stop to avoid fraying
- Single patient use recommended
- Can be wiped off or laundered

The Smart Glide rolls are delivered with a special knife, the Smart Glide Cutter, which makes it easy cutting out the exact size of slide sheet needed for the actual care task thus avoiding waste.



Smart Glide rolls

284472 M : width 100 cm, rolls of 100 m

284473 L : width 140 cm, rolls of 100 m

Smart Glide Precut

561500 : 5 boxes of 20 pcs, 100 x 125 cm

561500-1 : 1 box of 20 pcs, 100 x 125 cm

561501 : 5 boxes of 20 pcs, 100 x 145 cm

561501-1 : 1 box of 20 pcs, 100 x 145 cm

| Time to care |

V. Guldmann A/S
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com