



Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.guldmann.net

Cross hanger, 880 lbs Technical specifications

Functions

Lifting capacity, max. : 880 lbs (400 kg)
Suspension points : 4

Weight and Materials

Weight, totally : 10 lbs (4,6 kg)

Lifting hanger : Rustproof steel, powder coated

Measurements

A : 7 ³/₄" (192 mm)
B : 12" (309 mm)
C : 20 ¹/₂" (520 mm)

Accessories/Miscellaneous

Lifting sling : See brochure

Marking


The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device Class I.

Item Number/Model

561610 : Cross hanger, 880 lbs

Classification

according to ISO 9999 : 12 36 03



Time to care

Guldmann™