



Combi-lock


– Easy and safe to use

The Combi-Lock is used to connect and lock the rails together when two rail systems are to be combined. It combines the rail system between rooms, e.g. bedroom and bathroom to achieve an easy and comfortable transfer of a person.

The Combi-Lock is easy and safe to use. When the traverse rail is positioned correctly in front of the fixed rail, it locks - and the rails are automatically connected.

Connecting rail systems with the Combi-Lock makes it easy to design a complete rail system that meets any need.

The sleek design of the Combi-Lock blends easily with the surrounding architecture. It is easy to clean and maintain.



Easy and safe to use

- Intuitive and simple operating
- Effortless to use
- Automatic connection when the rails are correctly positioned
- High degree of safety

Nice and neat design

- All mechanical parts encased
- No hand control necessary
- Standard SWL 825 lbs

Improved Installation

- No need for transformer
- No need for power hook up / preparation
- Fast installation

Patent pending.



| Time to care |

V. Guldmann A/S
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com

Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.guldmann.net

**Guldmann Care-Lift
Solutions ULC**
Canada
Toll Free 1 877 339 7779
info@guldmann.net
www.guldmann.com