



GL6 Advanced

Mobile lifter

**Safer transfers.
More time to care.**

The GL6 is a versatile mobile lifter designed for hospitals, care homes and rehabilitation settings. It supports safe and efficient transfers in everyday care – from beds and chairs to wheelchairs and the floor.

GL6 combines dependable performance with thoughtful design – helping caregivers perform transfers with confidence and spend more time on care.

Designed for everyday care

GL6 combines ease of use with the flexibility needed in busy healthcare environments.

- Ergonomic soft-grip handle for comfortable and secure handling
- Flexible quick-release system allows the lifting hanger to move freely in all directions for easy positioning
- Electric leg-width adjustment for easy access around beds, chairs and wheelchairs
- Low-friction castors for smooth manoeuvring

- Compact turning radius makes it easier to move around furniture and in confined spaces
- Low leg profile helps position the lifter close to beds and other equipment
- Wide leg opening provides access around wheelchairs and other equipment
- Available with 155 kg or 205 kg lifting capacity



GL6 Advanced

For care environments where faster lifting and long battery life are important

- Lifting capacity: 155 kg or 205 kg
- Lifting speed: 36 mm/sec
- Long-life lithium-ion batteries
- On-board or off-board charging
- Electric leg-width adjustment
- Manual lowering as standard
- Low friction castors
– Ø75 mm or Ø100 mm

Flexible across care situations

GL6 provides a floor-based lifting solution when a ceiling hoist is not available or compatible. It can be used for both planned transfers and unforeseen lifting situations, including lifting from the floor.

- Transfers between floor, bed, chair and wheelchair
- Rehabilitation and treatment activities
- Low and high reach for a wide range of transfer situations
- Excellent reach across the bed
- Horizontal lifting with Guldmann stretcher-style accessories

Safety when it matters

GL6 is designed with accessible safety functions and clear, intuitive controls to support confident operation.

- Manual lowering as standard enables controlled lowering if power is unavailable
- Emergency stop immediately stops powered movement when activated
- Clear indicators and simple controls support straightforward operation
- Battery on/off function helps ensure the lifter is ready when needed and switched off when not in use

| Time to care |

V. Guldmann A/S
Tel. +45 8741 3100
info@guldmann.com
www.guldmann.com

